



# Lumit™ Immunoassay Cellular System Application Note


## Cellular Pathway Analysis Series

### Phospho 4E-BP1 (Ser 65)

#### Lumit™ Immunoassay Cellular System:

The Lumit™ Immunoassay Cellular System is a homogeneous bioluminescent assay that measures levels of target proteins in cell lysates when used with the appropriate primary antibody pairs (1). It combines immunodetection and NanoLuc Binary Technology (NanoBiT®) (2). In the Lumit™ Immunoassay Cellular System, NanoBiT® subunits (SmBiT and LgBiT) are conjugated to a pair of secondary antibodies against two different species (anti-rabbit, anti-mouse, or anti-goat). Seeded cells are lysed in multi-well plates using a Lumit™ compatible lysis solution and the target protein is detected by adding an antibody mix containing two primary antibodies against the target protein along with Lumit™ secondary antibodies. Binding of the primary/Lumit™ secondary antibody complexes to their corresponding epitopes brings NanoBiT® subunits into proximity to form an active NanoLuc® luciferase that makes light in proportion to the amount of the target protein (Fig. 1).


1. Hwang, B. *et al.* (2020) A homogeneous bioluminescent immunoassay approach to probe cellular signaling pathway regulation. *Commun Biol* 3, 8. doi:10.1038/s42003-019-0723-9.
2. Dixon, A. S. *et al.* (2016) NanoLuc Complementation Reporter Optimized for Accurate Measurement of Protein Interactions in Cells. *ACS Chem Biol* 11, 400-408.



**Figure 1. Illustration of Lumit™ Cellular Immunoassay.** When the primary antibody pair includes a phospho-specific antibody, the luminescence reflects the level of the target protein phosphorylation (top panel). To detect total protein level, the same concept is used except both primary antibodies recognize non-phosphorylated epitopes on the protein (bottom panel). The luminescent signal generated is measured using a luminometer.

#### Phospho 4E-BP1 (Ser 65) Immunoassay:

Upon activation of mTOR pathway with insulin, 4E-BP1 is phosphorylated (Fig. 2). After lysis of the cell membrane, phospho 4E-BP1 (Ser 65) can be detected using the reagents in Lumit™ Immunoassay Cellular System – Set 2 in combination with the anti 4E-BP1 antibodies described in Table 1.




**Figure 2. Detection of phosphorylated 4E-BP1 using the Lumit™ Immunoassay Cellular System – Set 2.** 50,000 seeded MCF-7 cells were starved overnight. The cells were then left untreated or pretreated with Rapamycin compound (10nM, 1hr) before they were untreated or treated with Insulin (0.1µM) for 90min. Phospho 4E-BP1 levels were measured following Promega Technical Manual TM613 and using the primary antibody conditions described in Table 1.




# Lumit™ Immunoassay Cellular System Application Note

## Cellular Pathway Analysis Series

### A Activation of 4E-BP1 phosphorylation with Insulin



### B Inhibition of 4E-BP1 phosphorylation with Rapamycin



**Figure 3. Activation and Deactivation of mTOR pathway.** (A) 50,000 seeded MCF-7 cells were starved overnight. Then they were untreated or treated with various concentrations of Insulin for 90min before phospho 4E-BP1 was measured by Lumit™ Immunoassay Cellular System – Set 2 to determine the Insulin EC<sub>50</sub>. (B) After starvation, 50,000 seeded MCF-7 were pretreated with various concentrations of Rapamycin for 1hr and then treated with Insulin (0.1μM, 90min) before phospho 4E-BP1 was measured by Lumit™ Immunoassay Cellular System – Set 2 to determine the potency of the inhibitor (IC<sub>50</sub>).

### Lumit™ Immunoassay Cellular System Short Protocol

1. Add 10μl lysis solution to 40μl cells.
2. Incubate for 20min with shaking.
3. Add 50μl Antibody mix.
4. Incubate for 60-90 min.
5. Add 25μl of Lumit™ detection reagent.
6. Shake plate for 2min.
7. Read luminescence.

This is a quick reference protocol. For more details regarding cells and reagent preparation and detailed protocols see Lumit™ Immunoassay Cellular System Technical Manual TM613 at [www.promega.com/protocols](http://www.promega.com/protocols).

**Table 1.**

Antibody*	Target	Supplier	Cat. #	Working stock (μg/ml)
p-4E-BP1 (Mouse)	Ser65	Thermo Fisher Scientific	MA5-31843	50
4E-BP1 (Rabbit)	Total	Cell Signaling Technology	9644	50

\*Antibodies from other suppliers may work as well. They may need optimization following Promega Technical Manual TM613.



### Ordering Information:

Products	Size	Promega Cat. #
Lumit™ Immunoassay Cellular System – Set 2	100 assays	W1331
	1,000 assays	W1332
	10,000 assays	W1333