

Total B-cell lymphoma 6 protein (BCL-6)

Lumit™ Immunoassay Cellular System:

The Lumit™ Immunoassay Cellular System is a homogeneous bioluminescent assay that measures levels of target proteins in cell lysates when used with the appropriate primary antibody pairs (1). It combines immunodetection and NanoLuc Binary Technology (NanoBiT®) (2). In the Lumit™ Immunoassay Cellular System, NanoBiT® subunits (SmBiT and LgBiT) are conjugated to a pair of secondary antibodies against two different species (anti-rabbit, anti-mouse, or anti-goat). Seeded cells are lysed in multi-well plates using a Lumit™ compatible lysis solution and the target protein is detected by adding an antibody mix containing two primary antibodies against the target protein along with Lumit™ secondary antibodies. Binding of the primary/Lumit™ secondary antibody complexes to their corresponding epitopes brings NanoBiT® subunits into proximity to form an active NanoLuc® luciferase that makes light in proportion to the amount of the target protein (Fig. 1).

- Hwang, B. *et al.* (2020) A homogeneous bioluminescent immunoassay approach to probe cellular signaling pathway regulation. *Commun Biol* 3, 8. doi:10.1038/s42003-019-0723-9.
- Dixon, A. S. *et al.* (2016) NanoLuc Complementation Reporter Optimized for Accurate Measurement of Protein Interactions in Cells. *ACS Chem Biol* 11, 400-408.




Figure 1. Illustration of Lumit™ Cellular Immunoassay. When the primary antibody pair includes a phospho-specific antibody, the luminescence reflects the level of the target protein phosphorylation (top panel). To detect total protein level, the same concept is used except both primary antibodies recognize non-phosphorylated epitopes on the protein (bottom panel). The luminescent signal generated is measured using a luminometer.

Total BCL-6 Immunoassay:

Upon treatment of cells with small molecule degrader BI-3802, BCL-6 is degraded (Fig. 2). After lysis of the cell membrane, total BCL-6 can be detected using the reagents in Lumit™ Immunoassay Cellular System – Set 1 in combination with the anti-BCL-6 antibodies described in Table 1.





Figure 2. Detection of total BCL-6 protein using the Lumit™ Immunoassay Cellular System – Set 1. 100,000 seeded Daudi cells were untreated or treated with a BCL-6 degrader compound, BI-3802 or a non-degrader inhibitor BI-3812 (100nM) for 2 hours. Total BCL-6 (Panel A) or a non-related control target, IκBα (Panel B), levels were measured following Technical manual TM613 and using the primary antibody conditions described in Table 1.




Figure 3. Targeted degradation of BCL-6 with small molecule degrader. 100,000 seeded Daudi cells were treated with various concentrations of BI-3802 or BI-3812 for 2 hours before BCL-6 and IκBα were measured by Lumit™ Immunoassay Cellular System – Set 1 to determine the potency of the small molecules (DC₅₀).

Lumit™ Immunoassay Cellular System Short Protocol

1. Add 10μl lysis solution to 40μl cells.
2. Incubate for 20min with shaking.
3. Add 50μl Antibody mix.
4. Incubate for 60-90 min.
5. Add 25μl of Lumit™ detection reagent.
6. Shake plate for 2min.
7. Read luminescence.

This is a quick reference protocol. For more details regarding cells and reagent preparation and detailed protocols see Lumit™ Immunoassay Cellular System Technical Manual TM613 at www.promega.com/protocols.

Table 1.

Antibody*	Target	Supplier	Cat. #	Working stock (μg/ml)
BCL-6 (Rabbit)	Total	Abcam	ab33901	50
BCL-6 (Mouse)	Total	Abcam	ab238044	50

*Antibodies from other suppliers may work as well. They may need optimization following Promega Technical Manual TM613.



Ordering Information:

Products	Size	Promega Cat. #
Lumit™ Immunoassay Cellular System – Set 1	100 assays	W1201
	1,000 assays	W1202
	10,000 assays	W1203