



Lumit™ Immunoassay Cellular System Application Note Cellular Pathway Analysis Series

Total Estrogen Receptor (ER)

Lumit™ Immunoassay Cellular System:

The Lumit™ Immunoassay Cellular System is a homogeneous bioluminescent assay that measures levels of target proteins in cell lysates when used with the appropriate primary antibody pairs (1). It combines immunodetection and NanoLuc Binary Technology (NanoBiT®) (2). In the Lumit™ Immunoassay Cellular System, NanoBiT® subunits (SmBiT and LgBiT) are conjugated to a pair of secondary antibodies against two different species (anti-rabbit, anti-mouse, or anti-goat). Seeded cells are lysed in multi-well plates using a Lumit™ compatible lysis solution and the target protein is detected by adding an antibody mix containing two primary antibodies against the target protein along with Lumit™ secondary antibodies. Binding of the primary/Lumit™ secondary antibody complexes to their corresponding epitopes brings NanoBiT® subunits into proximity to form an active NanoLuc® luciferase that makes light in proportion to the amount of the target protein (Fig. 1).

1. Hwang, B. *et al.* (2020) A homogeneous bioluminescent immunoassay approach to probe cellular signaling pathway regulation. *Commun Biol* 3, 8. doi:10.1038/s42003-019-0723-9.
2. Dixon, A. S. *et al.* (2016) NanoLuc Complementation Reporter Optimized for Accurate Measurement of Protein Interactions in Cells. *ACS Chem Biol* 11, 400-408.




Figure 1. Illustration of Lumit™ Cellular Immunoassay. When the primary antibody pair includes a phospho-specific antibody, the luminescence reflects the level of the target protein phosphorylation (top panel). To detect total protein level, the same concept is used except both primary antibodies recognize non-phosphorylated epitopes on the protein (bottom panel). The luminescent signal generated is measured using a luminometer.

Total Estrogen Receptor (ER) Immunoassay:

Upon treatment of cells with small molecule degrader Fulvestrant, ER is degraded (Fig. 2). After lysis of the cell membrane, total ER can be detected using the reagents in **Lumit™ Immunoassay Cellular System – Set 1** in combination with the anti ER antibodies described in Table 1.




Figure 2. Detection of total ER protein using the Lumit™ Immunoassay Cellular System – Set 1. 50,000 seeded MCF-7 cells were starved overnight. The cells were then untreated or treated with Fulvestrant or Tamoxifen (50nM) for 4 hours. Total Estrogen receptor levels were measured following Promega Technical Manual TM613 and using the primary antibody conditions described in Table 1.

Activation of Estrogen Receptor Degradation with Fulvestrant




Figure 3. Targeted degradation of Estrogen Receptor with small molecule degrader. After starvation overnight, 50,000 seeded MCF-7 cells were treated with various concentrations of Fulvestrant or Tamoxifen for 4 hours before ER (A) or ER and IκBα (B) were measured by Lumit™ Immunoassay Cellular System – Set 1 to determine the potency of the small molecules (DC₅₀).

Lumit™ Immunoassay Cellular System Short Protocol

1. Add 10μl lysis solution to 40μl cells.
2. Incubate for 20min with shaking.
3. Add 50μl Antibody mix.
4. Incubate for 60-90 min.
5. Add 25μl of Lumit™ detection reagent.
6. Shake plate for 2min.
7. Read luminescence.

This is a quick reference protocol. For more details regarding cells and reagent preparation and detailed protocols see Lumit™ Immunoassay Cellular System Technical Manual TM613 at www.promega.com/protocols.

Table 1.

Antibody*	Target	Supplier	Cat. #	Working stock (μg/ml)
ER (Rabbit)	Total	Thermo Fisher Scientific	MA5-14501	Supplied at stock concentration
ER (Mouse)	Total	Thermo Fisher Scientific	MA5-13191	50

*Antibodies from other suppliers may work as well. They may need optimization following Promega Technical Manual TM613.



Ordering Information:

Products	Size	Promega Cat. #
Lumit™ Immunoassay Cellular System – Set 1	100 assays	W1201
	1,000 assays	W1202
	10,000 assays	W1203