

SARS-CoV-2 3CLpro and PLpro Luminescent Assays Application Note

Luminescent Protease Assays:

Luminescent protease assays employ peptide aminoluciferin (peptide-aLuc) substrates in a homogeneous format that is ideally configured in multi-well plates and used for screening and characterizing protease inhibitors. In a 1st reaction a protease cleaves a substrate peptide moiety to release aminoluciferin (aLuc) that accumulates and drives a 2nd reaction with a luciferase that produces light in proportion to protease activity (Fig.1). The system uses a highly stabilized luciferase (UltraGlo™ Luciferase) in a Luciferin Detection Reagent that produces glow-style luminescence with a typical $t_{1/2} \geq 2$ hours.






Figure 1. Luminescent protease assay scheme. Protease and test article or vehicle are combined in an opaque white multi-well plate and a reaction is initiated by addition of a peptide-aLuc substrate. The reaction is stopped, and luminescence is initiated by adding Luciferin Detection Reagent. Signal is recorded on a plate-reading luminometer. Inhibitors are identified as test articles that reduce light output.

SARS CoV-2 3CLpro and PLpro luminescent assays:

3CLpro (a.k.a. Main Protease or Mpro) is a chymotrypsin-like protease and PLpro a papain-like protease. Both are encoded in the SARS-CoV-2 genome and play essential roles in the lifecycle of this virus. With the aLuc moiety shown in blue, the luminogenic 3CLPro Substrate is Ac-TSTKLQ-aLuc:



and the luminogenic PLpro Substrate is Z-RLRGG-aLuc:



Figures 2-4 show data from reactions of either recombinant SARS-CoV-2 3CLpro with Ac-TSTKLQ-aLuc or recombinant SARS-CoV-2 PLpro with Z-RLRGG-aLuc.




Figure 2. SARS-CoV-2 3CLpro enzyme assay with Luminogenic Substrate. 3CLpro (BPS, Cat# 100823) was combined with 20μM Ac-TSTKLQ-aLuc, 50mM HEPES (pH 7.2), 10mM DTT, and 0.1mM EDTA in a 50μl reaction volume in an opaque white 96-well plate and incubated for 1 hour at 37°C. Reactions were then terminated by adding 50μL of Luciferin Detection Reagent (Promega Cat.# V8920/V8921) and after 20 minutes at room temperature (20° – 25°C) luminescence was recorded on a GloMax® luminometer (Promega Cat.# GM2000). See step-by-step test compound screening protocol at the end of this document.




Figure 3. Detection of SARS-CoV-2 3CLpro inhibition with luminescent assay. 0.2μg of untagged* SARS-CoV-2 3CLpro (BPS Bioscience cat# 100823) was combined with 20μM Ac-TSTKLQ-aLuc, the 3CLpro inhibitor GC376 (AmBeed Cat.# A723661, diluted from 4mM stock DMSO solution; see Gurard-Levin et al Antiviral Research 182 (2020) 104924), 50mM HEPES (pH 7.2), 10mM DTT, and 0.1mM EDTA in a 50μl reaction volume in an opaque white 96-well plate, and incubated for 60 minutes at 37°C. Reactions were then terminated by adding 50μL of Luciferin Detection Reagent (Promega Cat.# V8920/V8921) and after 20 minutes at room temperature (20° – 25°C) luminescence was recorded on a GloMax® luminometer (Promega Cat.# GM2000). See step-by-step test compound screening protocol at the end of this document.

*Note: A 6His-tagged preparation of purified recombinant SARS-CoV-2 3CLpro that we tested displayed <1% of the activity of this untagged version and was insensitive to various reported 3CLpro inhibitors. This was consistent with a previous report that affinity tags depressed the activity of the highly homologous SARS-CoV 3CLpro (Grum-Tokers et al, Virus Res 133 (2008) 63–73).




Figure 4. SARS-CoV-2 PLpro enzyme assay with Luminogenic Substrate. GST-PLpro (R&D Systems, Cat# E-611-050) was combined with 20 μ M Z-RLRGG-aLuc, 50mM HEPES (pH 7.2), 10mM DTT, and 0.1mM EDTA in a 50 μ l reaction volume in an opaque white 96-well plate and incubated for 30 minutes at 20° – 25°C. Reactions were then terminated by adding 50 μ L of Luciferin Detection Reagent (Promega Cat.# V8920V8921) and after 20 minutes at room temperature (20° – 25°C) luminescence was recorded on a GloMax® luminometer (Promega Cat.# GM2000). See step-by-step test compound screening protocol at the end of this document.




Figure 5. Detection of SARS-CoV-2 PLpro inhibition with luminescent assay. 0.5nM SARS-CoV-2 GST-PLpro (R&D Systems, Cat# E-611-050) was combined with 20 μ M Z-RLRGG-aLuc, the PLpro inhibitor HY-17542 (MedChemExpress, diluted from 100mM stock DMSO solution; Loffredo et al <https://doi.org/10.1101/2020.07.15.203059>), 50mM HEPES (pH 7.2), 10mM DTT, and 0.1mM EDTA in a 50 μ l reaction volume in an opaque white 96-well plate, and incubated for 30 minutes at 20° – 25°C. Reactions were then terminated by adding 50 μ L of Luciferin Detection Reagent (Promega Cat.# V8920/V8921) and after 20 minutes at room temperature luminescence was recorded on a GloMax® luminometer (Promega Cat.# GM2000). See step-by-step test compound screening protocol at the end of this document.

The luminogenic SARS-CoV-2 protease substrates are currently available as early access materials from Promega along with the Promega catalogue items described in this Application Note as follows:

Product	Size	Cat. #
Luciferin Detection Reagent	10ml	V8920
Luciferin Detection Reagent	50ml	V8921
Luminogenic PLpro Substrate	please enquire*	
Luminogenic 3CLpro Substrate	please enquire*	
GloMax® Navigator Microplate Luminometer		GM2000

*Michael Curtin, Ph.D., Product Manager

Email: michael.curtin@promega.com

Phone: 608 298 4650

Step-by-step compound screening protocols

SARS-CoV-2 PLpro Assay

Materials

- The PLpro substrate is supplied at 4mM in 0.5M HEPES, pH 7.2 as a custom material from Promega Corp. Please inquire:
Michael Curtin, Ph.D., Product Manager
Email: michael.curtin@promega.com
Phone: 608 298 4650
- Assay buffer: 50mM HEPES pH 7.2, 10mM DTT, and 0.1mM EDTA (prepared by user)
- Luciferin Detection Reagent (Promega Cat.# V8920 or V8921)
- Recombinant PLpro sources:
 - R&D Systems Cat.# E-611-050
 - AcroBiosystems Cat.# PAE-C5148
- Opaque white 96 well plates (e.g., Corning Cat.# 3912)
- Plate reading luminometer (e.g., Promega GloMax, Cat.# GM2000)

Preparing reagents:

- Prepare 2X PLpro substrate solution at RT^o: Dilute PLpro substrate to 40μM in assay buffer.
- Prepare 4X PLpro enzyme solution on ice: Dilute PLpro recombinant enzyme to 2nM in assay buffer (see Figure 3 to consider enzyme concentration adjustments).
- Prepare 4X test compound solutions in assay buffer at RT^o.
- Prepare Luciferin Detection Reagent at RT^o: this reagent is supplied in 2 components, a lyophilized preparation, and a reconstitution buffer.
 - Add the entire contents of the reconstitution buffer to the Luciferin Detection Reagent lyophilized cake. Mix thoroughly but gentle to avoid forming bubbles (store unused portion at -20°C).

Performing the assay:

1. Dispense 25μL 2X PLpro substrate solution into an opaque white 96-well plate.
2. Add 12.5μL 4X test compound solutions to substrate solution in the 96-well plate.

3. Initiate reactions by adding 12.5 μ L 4X PLpro solution to substrate solution in the 96 well plate.
4. Mix plate and incubate for 30 minutes at room temperature (20°-25°C)
5. Add 50 μ L Luciferin Detection Reagent to each well to stop reactions and initiate luminescent signals. Allow 10 minutes for signal stabilization.
6. Read luminescence on a plate reading luminometer (e.g. Promega GloMax, Cat.# GM2000)

SARS-CoV-2 3CLpro Assay

Materials

- The 3CLpro substrate is supplied at 4mM in 0.5M HEPES, pH 7.2 as a custom material from Promega Corp. Please inquire:
Michael Curtin, Ph.D., Product Manager
Email: michael.curtin@promega.com
Phone: 608 298 4650
- Assay buffer: 50mM HEPES pH 7.2, 10mM DTT, and 0.1mM EDTA (prepared by user)
- Luciferin Detection Reagent: Promega Cat.# V8920 or V8921
- Recombinant 3CLpro: BPS Bioscience Cat.# 100823
- Opaque white 96 well plates (e.g., Corning, Cat.# 3912)
- Plate reading luminometer (e.g., Promega GloMax, Cat.# GM2000)

Preparing reagents:

- Prepare a 2X 3CLpro substrate solution (40 μ M) at RT°: Dilute 4mM 3CLpro substrate to 40 μ M in assay buffer.
- Prepare 4X 3CLpro enzyme solution on ice: Dilute 3CLpro untagged recombinant enzyme to 16 μ g/mL in assay buffer (see Figure 2 to consider enzyme concentration adjustments).
- Prepare 4X test compound solutions in assay buffer at RT°.
- Prepare Luciferin Detection Reagent at RT°: this reagent is supplied in 2 components, a lyophilized preparation, and a reconstitution buffer.
 - Add the entire contents of the reconstitution buffer to the Luciferin Detection Reagent lyophilized cake. Mix thoroughly but gentle to avoid forming bubbles (store unused portion at -20°C).

Performing the assay:

1. Dispense 25 μ L 2X 3CLpro substrate solution into an opaque white 96-well plate.
2. Add 12.5 μ L 4X test compound solutions to substrate solution in the 96-well plate.
3. Initiate reactions by adding 12.5 μ L 4X 3CLpro enzyme solution to substrate solution in the 96 well plate.
4. Mix plate and incubate for 1 hour at 37°C.
5. Add 50 μ L Luciferin Detection Reagent to each well to stop reactions and initiate luminescent signals. Allow 20 minutes for signal stabilization.
6. Read luminescence on a plate reading luminometer (e.g., Promega GloMax, Cat.# GM2000).