

Purification of Viral RNA from Sputum with the ReliaPrep™ Viral Total Nucleic Acid Purification Kit, Custom

Purify viral RNA from sputum using the ReliaPrep™ Viral Total Nucleic Acid Purification Kit, Custom.

Kit: ReliaPrep™ Viral Total Nucleic Acid Purification Kit, Custom (Cat.# AX4820)

Analyses: RT-qPCR for detection of Respiratory Syncytial Virus (RSV) and Influenza B

Sample Type(s): Sputum

Input: 200µl

Materials Required:

- DTT, Molecular Grade (Cat.# V3151)
- Nuclease-Free Water (Cat.# P1193)
- PBS, pH 7.2 (Gibco Cat.# 20012027) or similar
- ReliaPrep™ Viral Total Nucleic Acid Purification Kit, Custom (Cat.# AX4820)
- Microcentrifuge
- 1.5ml tubes
- Heat block set to 56°C
- 100% Isopropanol

This protocol was developed by Promega Applications Scientists and is intended for research use only.

Users are responsible for determining suitability of the protocol for their application.

For further information, see Technical Manual TM330, available at:

www.promega.com/protocols

or contact Technical Services at: techserv@promega.com

Sputum Processing¹:

1. Weigh appropriate amount of DTT, and rehydrate in Nuclease-Free Water to 500mM final concentration. Mix gently by pipetting to dissolve. DTT must be freshly made.
2. Prepare a 1:51 dilution of the 500mM DTT in PBS, pH 7.2. For example, add 100µl of 500mM DTT to 5.0ml of PBS.
3. Add an equal volume of diluted DTT in PBS to the sputum sample.
4. Incubate at room temperature with intermittent mixing by inversion until liquified or for a maximum of 30 minutes.
5. Proceed with purification.


Purification Protocol:

Binding Buffer (BBA) provided with this kit is not required for processing sputum samples.

1. Add 20µl of Proteinase K (PK) Solution to each 1.5ml tube.
2. Transfer 200µl of DTT-treated sputum to each 1.5ml tube.
3. Add 200µl of Cell Lysis Buffer (CLD) to each tube. Vortex for 10 seconds.

4. Incubate samples at 56°C for 10 minutes.
5. While samples are incubating, place a ReliaPrep™ Binding Column into an empty Collection Tube.
6. Remove the tube from the heat block. Add 250µl of 100% Isopropanol. Vortex for 10 seconds.
7. Transfer the tube contents to the ReliaPrep™ Binding Column, cap it and place it in a microcentrifuge.
8. Centrifuge for 1 minute at maximum speed.
9. Remove the collection tube containing flowthrough, and discard the liquid.
10. Place the binding column into a new collection tube.
11. Add 500µl of Column Wash Solution (CWD) to the column, and centrifuge for 3 minutes at maximum speed. Discard the flowthrough.
12. Repeat Step 11 twice for a total of three washes.
13. Place the column in a clean 1.5ml tube.
14. Add 60µl of Nuclease-Free Water to the column. Centrifuge for 1 minute at maximum speed.
15. Save the eluate, and discard the ReliaPrep™ Binding Column.

Results:



Detection of RSV and Influenza B RNA extracted from sputum. Sputum was treated with diluted DTT in PBS for 30 minutes. RSV A and Influenza B (Hong Kong) virus were reconstituted from Helix Elite™ Inactivated Standard Inactivated Influenza A/B and Respiratory Syncytial Virus (Microbiologics Cat.# HE0044N) and spiked into treated sputum. High virus sample contains approximately 2×10^5 copies of Influenza B and RSV A per 200µl sample. Low virus sample is a 1:100 dilution of the high virus sample in treated sputum. 200µl of the spiked sputum was processed with the ReliaPrep™ Viral Total Nucleic Acid Purification Kit, Custom as described above. Following nucleic acid purification, presence of RSV A and Influenza B was detected by RT-qPCR using GoTaq® 1-Step Probe qPCR System (Cat.# A6121). Each reaction contained 5µl of eluate with 12.5µl of the GoTaq® Probe qPCR Master Mix with dUTP, 0.5µl of GoScript™ RT Mix for 1-Step RT-qPCR, 1000nM forward and reverse primers and 200nM probe for RSV² or Influenza B³, and Nuclease-Free Water added to a final volume of 25µl. 1-step RT-qPCR thermal cycling was as follows³: reverse transcription at 50°C for 30 minutes, hot-start activation at 95°C for 2 minutes, and then 45 cycles of denaturation at 95°C for 15 seconds and annealing/extension at 55°C for 30 seconds, with signal acquisition during the annealing/extension stage of cycling. Data represent the average of duplicate purifications amplified in duplicate. Error bars indicate the standard deviation. Influenza B was detected in one of four amplifications of RNA extracted from the low virus sample.

References:

1. Processing of Sputum Specimens for Nucleic Acid Extraction, Centers for Disease Control <https://www.cdc.gov/coronavirus/2019-ncov/downloads/processing-sputum-specimens.pdf> Accessed 3/12/2020.
2. Fry, A.M., *et al.*, (2010) The Burden of Hospitalized Lower Respiratory Tract Infection due to Respiratory Syncytial Virus in Rural Thailand, *PLoS One*. **5**, e15098.
3. Selvaraju, S.B., *et al.*, (2010) Evaluation of Three Influenza A and B Real-Time Reverse Transcription-PCR Assays and a New 2009 H1N1 Assay for Detection of Influenza Viruses, *Journal of Clinical Microbiology*. **48**, 3870-3875.