



## Technical Appendix

### Composition of Solutions

**Note:** Solution compositions for the Cloning Enzyme 10X reaction buffers are listed in Table 15.

#### Antibiotic Stock Solutions

**Ampicillin:** 100mg/ml in deionized water (filter-sterilized).

**Tetracycline:** 10mg/ml in 80% ethanol.

**Kanamycin:** 30mg/ml kanamycin sulfate in deionized water (filter-sterilized).

**Chloramphenicol:** 20mg/ml in 80% ethanol.

Store at -20°C. **Note:** Cell growth in liquid culture is inhibited by tetracycline concentrations greater than 10mg/ml.

#### Cell Lysis Solution

0.2N NaOH, 1% SDS in deionized water. Prepare fresh for each use.

#### Cell Resuspension Solution

25mM	Tris-HCl (pH 8.0)
10mM	EDTA
50mM	glucose

#### IPTG Stock Solution (0.1M)

1.2g	isopropyl β-D-thiogalactopyranoside (IPTG) (Cat.# V3951)
------	--

Add deionized water to 50ml final volume. Filter-sterilize (0.2µm) and store at 4°C.

#### LB Medium

10g	Bacto®-tryptone
5g	Bacto®-yeast extract
5g	NaCl

Add deionized water to approximately 1L. Adjust pH to 7.5 with 10N NaOH and autoclave.

#### LB/Antibiotic Plates (1L)

Add 15g of agar to 1 liter of LB medium and autoclave. Allow the medium to cool to 55°C before adding antibiotic to the specified final concentration (ampicillin: 100µg/ml; tetracycline: 12.5µg/ml; kanamycin: 30µg/ml; chloramphenicol: 20µg/ml). Pour 30–35ml of medium into 85mm petri dishes. If necessary, flame the surface of the medium with a Bunsen burner to eliminate bubbles. Let the agar harden overnight. Store at 4°C for <1 month. Tetracycline is light-sensitive; LB/tetracycline plates should be covered with foil.

#### M-9 Plates

6g	Na <sub>2</sub> HPO <sub>4</sub>
3g	KH <sub>2</sub> PO <sub>4</sub>
0.5g	NaCl
1g	NH <sub>4</sub> Cl
15g	agar

Add deionized water to approximately 1L. Adjust pH to 7.4 with 10N NaOH. Autoclave. Cool to 50°C. Then add:

2.0ml	1M MgSO <sub>4</sub>
0.1ml	1M CaCl <sub>2</sub>
10.0ml	20% glucose
1.0ml	1M thiamine-HCl

Filter the medium through a 0.2µm filter unit.

#### Phenol:Chloroform:Isoamyl Alcohol (25:24:1)

Mix equal parts of TE buffer and phenol and allow the phases to separate. Then mix 1 part of the lower, phenol phase with 1 part of chloroform:isoamyl alcohol (24:1).

#### Potassium Acetate Solution (pH 4.8)

Prepare 60ml of 5M potassium acetate. Add 11.5ml of glacial acetic acid and 28.5ml of deionized water. This solution will be 3M with respect to potassium and 5M with respect to acetate. Store at 4°C.

#### SOC Medium

2.0g	Bacto®-tryptone
0.5g	Bacto®-yeast extract
1ml	1M NaCl
0.25ml	1M KCl
1ml	2M Mg <sup>2+</sup> stock (1M MgCl <sub>2</sub> • 6H <sub>2</sub> O, 1M MgSO <sub>4</sub> • 7H <sub>2</sub> O), filter-sterilized
1ml	2M glucose, filter-sterilized

Add Bacto®-tryptone, Bacto®-yeast extract, NaCl and KCl to 97ml deionized water. Stir to dissolve. Autoclave and cool to room temperature. Add 2M Mg<sup>2+</sup> stock and 2M glucose stock, each to a final concentration of 20µM. Filter the complete medium through a 0.2µm filter unit. The pH should be 7.0.

#### TAE 10X Buffer

400mM	Tris base
200mM	Sodium acetate
10mM	EDTA

Adjust pH to 8.2 with glacial acetic acid.

#### TBE 10X Buffer

890mM	Tris base
890mM	Boric acid
19mM	EDTA

Adjust pH to 8.3.

#### TE Buffer

10mM	Tris-HCl (pH 8.0)
1mM	EDTA

#### X-Gal Stock Solution (50mg/ml)

100mg	5-bromo-4-chloro-3-indolyl β-D-galactoside (X-Gal)
-------	--

Dissolve in 2ml of N,N'-dimethylformamide. This stock solution is available from Promega (Cat.# V3941).

**Table 12. Cloning Enzymes: Promega's Quality Acceptance Criteria.**

Enzyme	Concentration (u/μl)	SDS-PAGE Purity	Contaminant Activity		
			Endonuclease: Supercoiled DNA, 1 μg	dsDNase: 50ng Radiolabeled DNA	RNase: 50ng Radiolabeled RNA
AgarACE® Enzyme (1)	0.15–0.30	>90%	Gel analysis, absence of visible nicking or cutting, 5u/4hr at 42°C	<1% release, 5u/4hr at 42°C	<3% release, 5u/4hr at 42°C
Alkaline Phosphatase, Calf Intestinal (2)	1	N/A	Gel analysis, absence of visible nicking or cutting, 5u/1hr at 37°C	<3% release, 5u/1hr at 37°C	<3% release, 5u/1hr at 37°C
RecA Protein (3,4)	N/A	>95%	Gel analysis, absence of visible nicking or cutting, 15μg/16hr at 37°C	<1% release, 5μg/4hr at 37°C	N/A
T4 DNA Ligase (2,5)	1–3	>90%	Gel analysis, absence of visible nicking or cutting, 5u/16hr at 37°C	<1% release, 20u/16hr at 37°C	<3% release, 20u/5hr at 37°C
T4 Polynucleotide Kinase (6)	5–10	>90%	Gel analysis, absence of visible nicking or cutting, 25u/5hr at 37°C	<3% release, 25u/3hr at 37°C	<3% release, 25u/3hr at 37°C
T4 RNA Ligase	9–12	>90%	20u/3hr at 37°C	<1% release, 20u/3hr at 37°C	<1% release, 20u/3hr at 37°C

N/A – Not Applicable

1. T-Vector, Blue/White, DNA recovery of 2μg of DNA ladder.

2. Blue/White Assay is performed.

3. ssDNA-dependent ATPase activity.

4. Strand exchange.

5. T-Vector, Lambda ligation and packaging, ssDNase.

6. Oligo 5' end-labeled to high specific activity.

**Table 13. Uses and Genotypes of Various Strains of *E. coli*.**

Strain	Genotype	Comments
BL21(DE3)	F <sup>-</sup> , ompT, hsdS <sub>B</sub> , (r <sub>B</sub> <sup>-</sup> , m <sub>B</sub> <sup>-</sup> ), dcm, gal, λ(DE3)	<i>In vivo</i> expression of T7 Promoter-driven genes (e.g., pET-5, pGEM®).
<b>BL21(DE3)pLysS<sup>(e)</sup></b>	F <sup>-</sup> , ompT, hsdS <sub>B</sub> , (r <sub>B</sub> <sup>-</sup> , m <sub>B</sub> <sup>-</sup> ), dcm, gal, λ(DE3), pLysS CmR	<i>In vivo</i> expression of T7 Promoter-driven genes (e.g., pET-5, pGEM®), pLysS provides tighter control of T7 RNA Pol expression.
<b>BMH 71-18 mutS*</b>	thi, supE, Δ(lac-proAB), [mutS::Tn10] [F <sup>-</sup> , proAB, laq1 <sup>α</sup> ZΔM15]	GeneEditor™ Site-Directed Mutagenesis System; mismatch repair deficient strain.
<b>CJ236</b>	F <sup>-</sup> , cat(pCJ105; M13 <sup>cm</sup> )dut, ung1, thi-1, relA1, spoT1, mcrA	Kunkel mutagenesis; dut(-), ung(-).
<b>C600</b>	thi-1, thr-1, leuB6, lacY1, tonA21, supE44	λgt10, Permissive host; allows both parental and recombinant phage to grow.
<b>C600hfl</b>	thi-1, thr-1, leuB6, lacY1, tonA21, supE44, hflA150, [chr::Tn10]	λgt10, Restrictive host; allows recombinant phage to grow in preference to parental phage.
DH1	recA1, endA1, gyrA96, thi-1, hsdR17 (r <sub>K</sub> <sup>-</sup> , m <sub>K</sub> <sup>+</sup> ), supE44, relA1	Parent of DH5. DH5 more efficiently transformed by large (40 - 60kb) plasmids.
DH5α <sup>®</sup>	θ80dlacZΔM15, recA1, endA1, gyrA96, thi-1, hsdR17 (r <sub>K</sub> <sup>-</sup> , m <sub>K</sub> <sup>+</sup> ), supE44, relA1, deoR, Δ(lacZYA-argF)U169	Common host for cDNA cloning; supports α-complementation, recA(-) and endA(-).
DH5α <sup>®</sup> F'	F', θ80dlacZΔM15, recA1, endA1, gyrA96, thi-1, hsdR17 (r <sub>K</sub> <sup>-</sup> , m <sub>K</sub> <sup>+</sup> ), supE44, relA1, deoR, Δ(lacZYA-argF)U169	Single-stranded DNA synthesis.
<b>ES1301 mutS*</b>	lacZ53, thyA36, rha-5, metB1, deoC, IN(rrnD-rrnE), [mutS201::Tn5]	Provided with the Altered Sites® Mutagenesis Systems; mismatch repair deficient (kan <sup>r</sup> ).
<b>HB101*</b>	thi-1, hsdS20 (r <sub>B</sub> <sup>-</sup> , m <sub>B</sub> <sup>-</sup> ), supE44, recA13, ara-14, leuB6, proA2, lacY1, rpsL20 (str <sup>r</sup> ), xyl-5, mtl-1, galK2	Common strain for propagating plasmids that do not allow α-complementation.
JM83	ara, Δ(lac-proAB), rpsL, θ80dlacZΔM15	Host for pUC plasmids; pBR322- recA(+), r <sub>K</sub> <sup>+</sup> , m <sub>K</sub> <sup>+</sup> .
JM101	supE, thi, Δ(lac-proAB), [F <sup>-</sup> , traD36, proAB, lacI <sup>α</sup> ZΔM15]	Host for M13mp vectors; recA(+), r <sub>K</sub> <sup>+</sup> .

Strains listed in boldface are available from Promega.

\* Indicates strains available as competent cells.

(continued)



**Table 13. Uses and Genotypes of Various Strains of *E. coli* (continued).**

Strain	Genotype	Comments
JM103	<i>endA1, hsdR, supE, sbcB15, thi-1, strA, Δ(lac-proAB)</i> , [F', <i>traD36, proAB, lacI<sup>q</sup>ZΔM15</i> ]	Host for M13mp vectors; <i>recA</i> (+), <i>r<sub>K</sub></i> <sup>+</sup> .
JM105	<i>endA1, thi, rpsL, sbcB15, hsdR4, Δ(lac-proAB)</i> , [F', <i>traD36, proAB, lacI<sup>q</sup>ZΔM15</i> ]	Host for M13mp vectors; <i>recA</i> (+), <i>r<sub>K</sub></i> <sup>+</sup> .
JM107	<i>endA1, thi, gyrA96, hsdR17 (r<sub>K</sub><sup>-</sup>, m<sub>K</sub><sup>+</sup>), relA1, supE44, Δ(lac-proAB)</i> , [F', <i>traD36, proAB, lacI<sup>q</sup>ZΔM15</i> ]	Host for M13mp vectors; <i>recA</i> (+), <i>r<sub>K</sub></i> <sup>+</sup> .
JM108	<i>endA1, recA1, gyrA96, thi, hsdR17 (r<sub>K</sub><sup>-</sup>, m<sub>K</sub><sup>+</sup>), relA1, supE44, Δ(lac-proAB)</i>	Recombination deficient mutation in <i>rec A</i> increases stability of plasmids.
<b>JM109*</b>	<i>endA1, recA1, gyrA96, thi, hsdR17 (r<sub>K</sub><sup>-</sup>, m<sub>K</sub><sup>+</sup>), relA1, supE44, Δ(lac-proAB)</i> , [F', <i>traD36, proAB, lacI<sup>q</sup>ZΔM15</i> ]	Common host for cloning. ssDNA synthesis; restriction(-), <i>rec</i> (-), allows $\alpha$ -complementation: Included with the majority of our plasmids; Maintain on M9 plates supplemented with thiamine (to maintain F' episome).
<b>JM109(DE3)</b>	<i>endA1, recA1, gyrA96, thi, hsdR17 (r<sub>K</sub><sup>-</sup>, m<sub>K</sub><sup>+</sup>), relA1, supE44, Δ(lac-proAB)</i> , [F', <i>traD36, proAB, lacI<sup>q</sup>ZΔM15</i> ], $\lambda$ (DE3)	<i>In vivo</i> expression of T7 Promoter-driven genes (e.g., pET-5, pGEM <sup>®</sup> ), allows $\alpha$ -complementation.
JM110	<i>rpsL, thr, leu, thi, hsdR17 (r<sub>K</sub><sup>-</sup>, m<sub>K</sub><sup>+</sup>), lacY, galK, galT, ara, tonA, tsx, dam, dcm, supE44, Δ(lac-proAB)</i> , [F', <i>traD36, proAB, lacI<sup>q</sup>ZΔM15</i> ]	<i>dam</i> (-), <i>dcm</i> (-) strain, allows $\alpha$ -complementation.
<b>LE392</b>	<i>hsdR514, (r<sub>K</sub><sup>-</sup>, m<sub>K</sub><sup>+</sup>), supE44, supF58, lacY1 or Δ(lacIZY)6, galK2, galT22, metB1, trpR55</i>	Genomic and cDNA cloning; restriction(-), <i>rec</i> (+), permissive host, no color selection, <i>lon</i> (+). Recommended (if no color selection needed) as primary strain for amplification of recombinant phage and screening of cDNA library with nucleic acid probe.
<b>KW251</b>	<i>supE44, galK2, galT22, metB1, hsdR2, mcrB1, mcrA, [argA81::Tn10], recD1014</i>	Genomic cloning; permissive host, used as alternative to LE392, <i>rec</i> (-) strain, restriction (-), <i>tet</i> <sup>r</sup> .
MB408	<i>recF, recB21, recC22, sbcB15, hflA, hflB, hsdR, (tet<sup>r</sup>)</i>	
MC1061	F-, <i>araD139, Δ(ara-leu)7696, galE15, galK16, Δ(lac)X74, rpsL(Str), hsdR2(r<sub>K</sub><sup>-</sup>, m<sub>K</sub><sup>+</sup>), mcrA, mcrB1</i>	
<b>NM522</b>	<i>supE, thi, Δ(lac-proAB), Δhsd5 (r-, m-)</i> , [F', <i>proAB, lacI<sup>q</sup>ZΔM15</i> ]	ssDNA synthesis; restriction (-), <i>rec</i> (+), F'. Grow on M9 plates to maintain F' episome.
<b>NM538</b>	<i>supF, hsdR (r<sub>K</sub><sup>-</sup>, m<sub>K</sub><sup>+</sup>), trpR, lacY</i>	Genomic cloning; permissive host.
<b>NM539</b>	<i>supF, hsdR (r<sub>K</sub><sup>-</sup>, m<sub>K</sub><sup>+</sup>), lacY, (P2)</i>	Restrictive host; used for <i>Spi</i> selection of recombinant phage.
P2392	LE392 (P2)	
RR1	<i>hsdS20, (r<sub>B</sub><sup>-</sup>, m<sub>B</sub><sup>+</sup>), supE44, ara-14, proA2, rpsL20 (str), lacY1, galK2, xyl-5, mtl-1, supE44</i>	<i>Rec A</i> (+) version of HBI01.
$\chi$ 1776	<i>tonA53, dapD8, minA1, glnV44, supE44, Δ(gal-uvrB)40, min82, rib-2, gyrA25, thyA142, oms-2, metC65, oms-1, (tte-1), Δ(bioH-asd)29, cyc82, cycA1, hsdR2</i>	Debilitated strain used in early work with recombinant DNA.
XL1-Blue	<i>recA1, endA1, gyrA96, thi, hsdR17(r<sub>K</sub><sup>-</sup>, m<sub>K</sub><sup>+</sup>), supE44, relA1, lac, [F', <i>proAB<sup>+</sup>, lacI<sup>q</sup>ZΔM15, ::Tn10(Tet<sup>r</sup>)</i>]</i>	Common host for cloning.
SURE Cells	<i>E14-(McrA-) Δ(mcrCB-hsdSMR-mrr) 171 endA1, supE44, thi-1, gyrA96, relA1, lac, recB, recJ, sbcC, umuC::Tn5(kan<sup>r</sup>), uvrC, [F', <i>proAB, lacI<sup>q</sup>ZΔM15 [Tn10(Tet<sup>r</sup>)</i>]</i>	Increased stability of DNA containing inverted repeats or Z-DNA (Stratagene); F' episome- able to make ssDNA from phagemid DNA.
Y1088	$\Delta$ ( <i>lacU169</i> ), <i>supE, supF, hsdR (r<sub>K</sub><sup>-</sup>, m<sub>K</sub><sup>+</sup>), metB, trpR, tonA21, [proC::Tn5] (pMC9)</i>	Host for amplification of cDNA libraries. pMC9 confers <i>amp<sup>r</sup></i> , <i>tet<sup>r</sup></i> .
<b>Y1089</b>	$\Delta$ ( <i>lacU169</i> ), <i>proA<sup>+</sup>, Δ(lon), araD139, strA, hflA150, [chr:Tn10(tet<sup>r</sup>)</i> , (pMC9)	cDNA cloning- <i>lon</i> (-), <i>hflA150</i> (enhances lysogeny), pMC9- maintain on <i>amp/tet</i> plates: Used primarily for generation of preparative amounts of recombinant fusion protein. pMC9 confers <i>amp<sup>r</sup></i> , <i>tet<sup>r</sup></i> .
<b>Y1090</b>	$\Delta$ ( <i>lacU169</i> ), <i>proA<sup>+</sup>, Δ(lon), araA139, strA, supF, [trpC22::Tn10 (tet<sup>r</sup>)</i> , (pMC9), <i>hsdR (r<sub>K</sub><sup>-</sup>, m<sub>K</sub><sup>+</sup>)</i>	cDNA cloning. Useful for screening expression cDNA libraries. pMC9 confers <i>amp<sup>r</sup></i> , <i>tet<sup>r</sup></i> .

Strains listed in boldface are available from Promega.

\* Indicates strains available as competent cells.



**Promega**

Precision Design... for Life

www.promega.com

**Table 14. Applications of Promega's Competent Cells.**

Applications	Cells	Genotype	>1 x 10 <sup>8</sup> cfu/μg	>1 x 10 <sup>7</sup> cfu/μg	>1 x 10 <sup>6</sup> cfu/μg
Cloning T-Vector Cloning Blue/White Screening	JM109	<i>endA1, recA1, gyrA96, thi hsdR17</i> ( <i>r<sub>K</sub><sup>-</sup>, m<sub>K</sub><sup>+</sup></i> ), <i>relA1, supE44 Δ(lac- proAB)</i> , [ <i>F' traD36, proAB</i> <i>lacI<sup>q</sup>ZΔM15</i> ].	✓	✓	
Cloning	HB101	<i>F', thi-1, hsdS20 (r<sub>B</sub><sup>-</sup>, m<sub>B</sub><sup>-</sup>), supE44, recA13, ara-14, leuB6, proA2, lacY1, galK2 rpsL20, (str<sup>r</sup>), xyl-5, mtl-1</i>	✓	✓	
Site-Directed Mutagenesis	BMH 71-18 <i>mutS</i>	<i>thi, supE, Δ(lac-proAB), [mutS::Tn10], [F' proAB, lacI<sup>q</sup>ZΔM15].</i> Mutagenesis		✓	
	ES1301 <i>mutS</i>	<i>lacZ53, mutS201::Tn5, thyA36, rha-5, metB1, decC1, IN(rrnD- rrnE).</i>		✓	
Protein Expression	BL21 (DE3)pLysS <sup>(e)</sup>	<i>F', ompT, hsdS<sub>B</sub>(r<sub>B</sub><sup>-</sup>, m<sub>B</sub><sup>+</sup>), dcm, gal, (DE3), pLysS, Cmr<sup>r</sup></i>		✓	✓

**Table 15. Cloning Enzyme 10X Reaction Buffer Formulations.**

Enzyme	10X Formulation
T4 DNA Ligase*	300mM Tris-HCl (pH 7.8), 100mM MgCl <sub>2</sub> , 100mM DTT, 10mM ATP
T4 RNA Ligase	500mM Tris-HCl (pH 7.8), 100mM MgCl <sub>2</sub> , 50mM DTT, 10mM ATP
T4 Polynucleotide Kinase (PNK)	700mM Tris-HCl (pH 7.6), 100mM MgCl <sub>2</sub> , 50mM DTT
Calf Intestinal Alkaline Phosphatase (CIAP)	500mM Tris-HCl (pH 9.3), 10mM MgCl <sub>2</sub> , 1mM ZnCl <sub>2</sub> , 10mM spermidine

\* pGEM<sup>®</sup>-T and pGEM<sup>®</sup>-T Easy Vector Systems (Cat.# A1360, A1380, A3600, A3610) provide a 2X Rapid Ligation Buffer for T4 DNA Ligase with a formulation of: 60mM Tris-HCl (pH 7.8), 20mM MgCl<sub>2</sub>, 20mM DTT, 2mM ATP, 10% polyethylene glycol (MW 8000).

The Rapid Ligation Buffer enables performance of ligation reactions in as little as 1 hour.



**General References**

1. Sambrook, J., Fritsch, E.F. and Maniatis, T. (1989) *Molecular Cloning: A Laboratory Manual*, Cold Spring Harbor Press, Cold Spring Harbor, NY.
2. Anderson, S. *et al.* (1980) *Nucl. Acids Res.* **8**, 1731.
3. Challberg, M.D. and Englund, P.T. (1980) *Meth. Enzymol.* **65**, 39.
4. Perbal, B. (1988) In: *A Practical Guide to Molecular Cloning*, 2nd ed., John Wiley and Sons, NY.
5. Burd, J.F. and Wells, R.D. (1974) *J. Biol. Chem.* **249**, 7094.
6. Kaufman, D.L. and Evans, G.A. (1990) *BioTechniques* **9**, 304.
7. Eckert, K.A. and Kunkel, T.A. (1992) In: *PCR, A Practical Approach*, McPherson, M.J. *et al.* eds., Oxford University Press, NY.
8. Hanahan, D. (1985) In: *DNA Cloning*, Volume 1, Glover, D. ed., IRL Press Ltd., London, United Kingdom.
9. Mezei, L.M. and Storts, D.R. (1994) In: *PCR Technology: Current Innovations*, Griffin, H.G. and Griffin, A.M. eds. CRC Press, Boca Raton, Florida.
10. Robles, J. and Doers, M. (1994) *Promega Notes* **45**, 19.
11. Clark, J.M. (1988) *Nucl. Acids Res.* **19**, 1156.
12. Lundberg, K.S. *et al.* (1991) *Gene* **108**, 1.
13. Kong, H. *et al.* (1993) *J. Biol. Chem.* **268**, 1965.
14. Kobs, G. (1996) *Promega Notes* **55**, 28.
15. Haff, L. and Mezei, L. (1989) *Amplifications* **1**, 8.

**Additional Literature**

**Technical Manuals and Product Information Sheets**

- #9PIM180 T4 DNA Ligase
- #9PIM182 Calf Intestinal Alkaline Phosphatase
- #9PIM181 T4 DNA Polymerase
- #9PIM220 DNA Polymerase I Large (Klenow) Fragment
- #TM042 pGEM®-T and pGEM®-T Easy Vector Systems

**Promega Notes Articles**

- PN032 pGEM®-Zf Vector update: Considerations for optimal single-stranded DNA production
- PN045 pGEM®-T Vector Systems troubleshooting guide
- PN051 pGEM®-T Vector: Technically Speaking
- PN054 pGEM®-T Vector: Cloning of modified blunt-ended DNA fragments
- PN071 Rapid Ligation for the pGEM®-T and pGEM®-T Easy Vector Systems  
Cloning Blunt End *Pfu* DNA Polymerase-Generated PCR Fragments into pGEM®-T Vector



## Glossary

**5'-RACE:** Rapid Amplification of cDNA 5'-Ends.

**Activated Calf Thymus DNA:** Nicked and gapped dsDNA prepared by treatment with DNase I; a substrate for many DNA polymerases.

**Alkaline Phosphatase (AP):** An enzyme that catalyzes the removal of a phosphate group from a substrate. This property is used in colorimetric and chemiluminescent detection reagents.

**Alpha Complementation:** The process by which a functional  $\beta$ -galactosidase (*lacZ*) gene is generated when the *lacZ*  $\alpha$ -peptide (N-terminus) complements the  $\omega$  fragment of *lacZ* (C-terminus). In general, a cloning vector contributes the *lacZ*  $\alpha$ -peptide and the host bacterial strain provides the  $\omega$  fragment. See also blue/white cloning.

**Blue/White Cloning:** A technique used to identify recombinant (positive) clones in cloning experiments. Specially constructed cloning vectors will produce  $\beta$ -galactosidase by  $\alpha$ -complementation when transformed into the appropriate host strain after exposure to the chemical IPTG. This will cause the cells to exhibit a blue color in the presence of the substrate X-Gal. Foreign DNA inserted into the cloning vector will interrupt the *lacZ* gene, preventing  $\alpha$ -complementation and producing white colonies.

**Blunt Ends:** DNA ends that do not possess short, single-stranded overhangs. See also sticky ends.

**Cloning:** The production of a large number of identical DNA molecules from a single DNA molecule through replication of the DNA within a cell. Cloning also refers to the production of large numbers of identical cells from a single cell.

**Cohesive Ends:** Two DNA ends, in the same or different molecules, which have short, single-stranded overhangs that are complementary to one another. Cohesive ends allow comparatively efficient ligation of DNA molecules. See also sticky ends.

**Competent Cells:** Bacterial cells that are able to take in exogenous DNA.

**Distributive:** Enzyme dissociates from template after a single nucleotide addition.

**DNA (Deoxyribonucleic Acid):** A polymeric molecule composed of deoxyribonucleotide units joined in a specific sequence through the formation of 3'→5' phosphodiester bonds.

**DNase (Deoxyribonuclease):** An enzyme that breaks down DNA into small fragments or deoxyribonucleotides.

**DNA Ligase:** An enzyme that joins two DNA molecules (or two ends of the same DNA molecule) by forming a phosphodiester bond between the 3'- and 5'-ends.

**dNTP:** Deoxyribonucleotide 5'-triphosphate.

**ds:** Double-stranded.

**DTT:** Dithiothreitol.

**End-Labeling:** The addition of a labeled group (radioactive or nonradioactive) to the 5'- or 3'-end of DNA or RNA. This is typically accomplished by using a kinase to label the 5'-end, or a DNA polymerase or terminal deoxynucleotidyl transferase to label the 3'-end.

**Endonuclease:** An enzyme that hydrolyzes phosphodiester bonds at internal locations within a DNA or RNA molecule.

**Exonuclease:** An enzyme that hydrolyzes phosphodiester bonds at the ends of DNA or RNA molecules, resulting in the stepwise removal of nucleotides.

**Gene:** A segment of DNA that encodes a polypeptide, protein or RNA molecule.

**Genotype:** The specific genes (which may or may not be expressed) that are present in an organism.

**Heat-Inactivation:** Destroying an enzyme's activity by heating to a high temperature (typically 60–70°C) for an extended length of time. Not all enzymes can be heat-inactivated.

**In vitro:** A reaction or experiment performed in the absence of living cells, typically using conditions that attempt to mimic those found within cells.

**In vivo:** A reaction or experiment performed in a living organism or cell.

**kb:** Kilobase or kilobase pairs.

**$k_{cat}$ :** Maximum number of substrate molecules converted to products per active site per unit time.

**kDa:** KiloDalton.

**$K_m$ :** The Michaelis constant; the concentration of substrate that an enzyme can convert to product at half its maximal rate.

**Labeling:** A process in which nucleic acids or proteins are tagged with a radioactive or nonradioactive marker.

**Ligase:** An enzyme that catalyzes DNA or RNA linkage, generally splitting off a pyrophosphate group from ATP concurrently.

**Multiple Cloning Site (MCS):** The region of a DNA vector that contains unique restriction enzyme recognition sites into which foreign DNA can be inserted; also called a polylinker.

**Nuclease:** An enzyme that degrades nucleic acids.

**Nucleotide:** A molecule composed of an organic base, sugar and phosphate group, which constitutes the "building blocks" of nucleic acids (DNA and RNA).

**Oligonucleotide (Oligo):** A short (typically <50 nucleotides), single-stranded DNA or RNA molecule.

**PCR:** Polymerase Chain Reaction.

**PEG:** Polyethylene glycol.

**Phosphatase:** An enzyme that removes a phosphate group from a protein, nucleic acid or other molecule.

**$PP_i$ :** Inorganic pyrophosphate.

**Primer:** An oligonucleotide or short single-stranded nucleic acid that acts as a starting point for the synthesis of nucleic acids from a template.

**Promoter:** DNA sequence for the initiation of RNA transcription by RNA polymerase.

**RNA (Ribonucleic Acid):** A polymeric molecule composed of ribonucleotide units joined in a specific sequence through the formation of 3'→5' phosphodiester bonds.

**RNase (Ribonuclease):** An enzyme that breaks down RNA into smaller RNA fragments or ribonucleotides.

**rNTP:** Ribonucleotide 5'-triphosphate.

**ss:** Single-stranded.

**Sticky Ends:** Two DNA ends, in the same or different molecules, that have short, single-stranded overhangs that are complementary to one another. Sticky ends allow comparatively efficient ligation of DNA molecules. See also Cohesive Ends.

**TAE:** Tris Acetate EDTA.

**TCA:** Trichloroacetic acid.


**Terminator:** DNA sequence for the termination of RNA transcription by RNA polymerase.

**Transformation:** The process during which a plasmid DNA is inserted into a bacterial cell.

**Turnover Rate ( $k_{cat}$ ):** Maximum number of substrate molecules converted to products per active site per unit time.

**Vector:** A DNA molecule that can replicate within a host cell and that allows the insertion of foreign DNA sequences. Vectors commonly used may be plasmids, phagemids or bacteriophage.





**Cover Art** adapted from Figure 3 p. 139 in Alberto I. Roca and Michael M. Cox, "RecA Protein: Structure, Function, and Role in Recombinatorial DNA Repair," in *Progress in Nucleic Acid Research*, Volume 56, pp. 129-224, edited by Waldo E. Cohn and Kiviva Moldave. Reprinted by permission. Copyright © 1997 by Academic Press.

AgarACE, Altered Sites, Erase-a-Base, *fmol*, OmniBase, Packagene, pGEM, Prime-a-Gene, RiboClone, TagTrack and Wizard are trademarks of Promega Corporation and are registered with the U.S. Patent and Trademark Office.

GeneEditor, MULTI-CORE, pBEST<sub>luc</sub>, PinPoint and pTARGET and SoftLink are trademarks of Promega Corporation.

Bacto is a registered trademark of Difco Laboratories, Detroit, Michigan. Coomassie is a registered trademark of Imperial Chemical Industries, Ltd. DH5α is a registered trademark of Life Technologies, Inc. Ficoll is a registered trademark of Amersham Pharmacia Biotech AB. GTG, NuSieve, SeaKern and SeaPlaque are registered trademarks of FMC Corporation. SpeedVac is a registered trademark of Savant Instruments. SURE is a registered trademark of Stratagene. SYBR is a registered trademark of Molecular Probes, Inc. Tween is a registered trademark of ICI Americas, Inc.

©U.S. Pat. No. 5,869,310 and other patents.

©U.S. Pat. No. 5,658,548 and other patents.

©Licensed under U.S. Pat. No 5,075,430.

©U.S. Pat. No. 4,766,072.

©**Usage Restrictions for the T7 Expression System**

The T7 expression system is based on technology developed at Brookhaven National Laboratory under contract with the U.S. Department of Energy and is the subject of patents and patent applications assigned to Brookhaven Science Associates, LLC (BSA). This technology, including bacteria, phage and plasmids that carry the gene for T7 RNA polymerase, is to be used for academic or nonprofit laboratory or licensed commercial research purposes only. By accepting or using the T7 expression technology you agree to be bound by the following conditions set forth by BSA. The initial purchaser may refuse to accept the conditions of this notice by returning this product and the enclosed materials to Promega unused.

**Academic and NonProfit Laboratories**

No materials that contain the cloned gene for T7 RNA polymerase may be distributed further to third parties outside of your laboratory unless the recipient receives a copy of this assurance notice and agrees to be bound by its terms. This limitation applies to Bacterial Strains JM109(DE3) and BL21(DE3)pLysS, to the following Promega products that include these bacterial strains (pGEMEX@-1 and pGEMEX@-2 Vectors, pET-5 and pET-9 Vectors, and all pET Vector Systems) and to any derivatives thereof.

**Commercial Laboratories**

A license is required for any commercial use of the T7 expression system, including use of the T7 system for research purposes or for production purposes by any commercial entity. Information about commercial licenses may be obtained from the Licensing Office, Brookhaven National Laboratory, Upton, NY 11973, Telephone: 516-344-3312, FAX: 516-344-3729.

©Licensed under one or both of U.S. Pat. No. 5,487,993 and European Pat. No 0 550 693.

©U.S. Pat. No. 5,780,270.

©Covered under U.S. Pat. No. 5,168,062 assigned to the University of Iowa Research Foundation.

©For research purposes only. Not for diagnostic or therapeutic use. For nonresearch uses of the portion of the vector encoding the biotinylation sequence, please contact Promega Corporation for licensing information.

©U.S. Pat. No. 5,552,302 and other patents. *Inhibitors of Angiogenin*, which comprises a segment of human PRI, is the subject of U.S. Pat. No. 4,966,964 and other patents assigned to the President and Fellows of Harvard College and exclusively licensed to Promega Corporation.

©Licensed from Boehringer Mannheim GmbH under U.S. Pat. No. 4,804,748.

©The method of recombinant expression of *Coleoptera* luciferase is covered by U.S. Pat. Nos. 5,583,024, 5,674,713 and 5,700,673.

**NORTH AMERICA**

**Promega Headquarters**  
2800 Woods Hollow Road  
Madison, Wisconsin 53711-5399  
Toll Free in USA (800) 356-9526  
Toll Free FAX in USA (800) 356-1970  
Phone (608) 274-4330  
FAX (608) 277-2516  
**Fisher Scientific Canada**  
Toll Free in Canada (800) 234-7437

**EUROPE**

▲ Austria +800-77663422  
▲ Belgium 0800-18098  
Czech Republic 02 301 8177  
Denmark 44 94 88 22  
Finland 358 9 7552550  
▲ France 0800 48 79 99  
▲ Germany +800-77663422  
Greece 1 6436138  
Hungary 1 251 0344  
Ireland 018224222  
▲ Italy 800 69 1818  
▲ Luxembourg +31(0)71 532 4244

▲ The Netherlands 0800-0221910  
Norway 22715090  
Poland 58 3414726  
Portugal 1 3613620  
Russia 095 135 4206  
Slovak Republic 4217 4341 0365  
Slovenia 0 64 351 510  
Spain 93 404 52 14  
Sweden/Iceland 0346 83050  
▲ Switzerland 01 830 70 37  
▲ United Kingdom 0800 378994

**LATIN AMERICA**

Argentina/Uruguay 1 4383 3000  
Brazil 22715090  
bioBRAS (*Belo Horizonte*) 31 2919877  
Chile 22472976  
Promicro (*São Paulo*) 118690699  
Colombia 1 6112287  
Ecuador 32582483  
Mexico 5 281 4718  
Venezuela 2 265 0891

**PACIFIC ASIA**

▲ Australia 1 800 225 123  
▲ Beijing 10 68498287  
Joint Venture SP 21 6483 5136  
SABC Hong Kong 379 432 2510  
India 11 684 6565  
Indonesia 21 489 1718  
▲ Japan 03 3869 7981  
Korea (02) 478 5911  
Malaysia 3 7048600  
New Zealand 93664784

Pakistan 92 51 518 469  
Singapore 775 7284  
Taiwan 0223825378  
Thailand 2 683 0114  
Bio-Active Co., Ltd. Diagnostic Biotech 2 295 2532  
**MIDDLE EAST/AFRICA**  
Egypt 2 245 1785  
Israel 8 9477077  
Lebanon 1 249836  
South Africa 21 981 1560  
Turkey 216 385 8321

▲ Indicates Promega Office



© 1999 Promega Corporation. All Rights Reserved.  
Prices and specifications subject to change without prior notice.



**Promega**

Precision Design... for Life

www.promega.com

Rev.041599.eps