

RSK3 Kinase Assay

By Juliano Alves, Laurie Engel, Said A. Goueli, and Hicham Zegzouti, Promega Corporation

Scientific Background:

RSK3 is a member of the RSK (ribosomal S6 kinase) family that encodes a 733-amino-acid protein with a unique N-terminal region containing a putative nuclear localization signal (1). RSK3 mRNA is widely expressed and is activated by growth factors, serum and phorbol ester. Upon stimulation, RSK3 translocates to the cell nucleus and phosphorylates nuclear proteins. RSK3 can bind to ERK1/2 and this association increases the duration of RSK3 activation (2).

1. Zhao, Y. et al: RSK3 encodes a novel pp90rsk isoform with a unique N-terminal sequence: growth factor-stimulated kinase function and nuclear translocation. *Mol Cell Biol.* 1995 Aug;15(8):4353-63.
2. Roux, P P. et al: Phosphorylation of p90 ribosomal S6 kinase (RSK) regulates extracellular signal-regulated kinase docking and RSK activity. *Mol Cell Biol.* 2003 Jul;23(14):4796-804.

ADP-Glo™ Kinase Assay

Description

ADP-Glo™ Kinase Assay is a luminescent kinase assay that measures ADP formed from a kinase reaction; ADP is converted into ATP, which is converted into light by Ultra-Glo™ Luciferase (Fig. 1). The luminescent signal positively correlates with ADP amount (Fig. 2) and kinase activity (Fig. 3A). The assay is well suited for measuring the effects chemical compounds have on the activity of a broad range of purified kinases—making it ideal for both primary screening as well as kinase selectivity profiling (Fig. 3B). The ADP-Glo™ Kinase Assay can be used to monitor the activity of virtually any ADP-generating enzyme (e.g., kinase or ATPase) using up to 1mM ATP.




Figure 1. Principle of the ADP-Glo™ Kinase Assay. The ATP remaining after completion of the kinase reaction is depleted prior to an ADP to ATP conversion step and quantitation of the newly synthesized ATP using luciferase/luciferin reaction.




Figure 2. Linearity of the ADP-Glo Kinase Assay. ATP-to-ADP conversion curve was prepared at 10µM ATP+ADP concentration range. This standard curve is used to calculate the amount of ADP formed in the kinase reaction. Z' factors were determined using 200 replicates of each of the % conversions shown.



ADP-Glo™ Kinase Assay Application Note Ser/Thr Kinase Series

The following is only a short protocol. For detailed protocols on conversion curves, kinase assays and inhibitor screening, see Kinase Enzyme Systems Protocol at: <http://www.promega.com/KESProtocol>

Short Protocol

- Dilute enzyme, substrate, ATP and inhibitors in 1x kinase reaction buffer.
- Add to the wells of 384 low volume plate:
 - ✓ 1 μ l of inhibitor or (5% DMSO)
 - ✓ 2 μ l of enzyme (defined from table 1)
 - ✓ 2 μ l of substrate/ATP mix
- Incubate at room temperature for indicated time (See Figure 3).
- Add 5 μ l of ADP-Glo™ Reagent.
- Incubate at room temperature for 40 minutes.
- Add 10 μ l of Kinase Detection Reagent.
- Incubate at room temperature for 30 minutes.
- Record luminescence (Integration time 0.5-1 second).

Table 1. Enzyme Titration. Data are shown as relative light units (RLU) that directly correlate to the amount of ADP produced. The correlation between the % of ATP converted to ADP and corresponding signal to background ratio is indicated for each kinase amount.

Enzyme, ng	100	50	25	12.50	6.25	3.13	1.56	0.78	0.39	0.20	0.10	0
Luminescence	216,865	202,140	212,488	199,656	181,383	127,946	90,087	47,511	26,593	16,286	9,241	1,701
S/B	127	119	125	117	107	75	53	28	16	10	5	1
% Conversion	77	72	76	71	64	45	31	16	8	4	2	0




Figure 3. RSK3 Kinase Assay Development. (A) RSK3 enzyme was titrated using 10 μ M ATP and the luminescence signal generated from each of the amounts of the enzyme is shown. (B) Inhibitor dose response was created using 1ng of RSK3 to determine the potency of the inhibitor (IC_{50}).

Ordering Information:



Products	Size	Cat. #
RSK3 Kinase Enzyme System	10 μ g	VA7549
	1mg	VA7550
ADP-Glo™ + RSK3 Kinase Enzyme System	1 Each	VA7551