

TECHNICAL MANUAL

Raji (LDH-HiBiT) TCK Bioassays

Instructions for Use of Products

JA1311, JA1315, JA1221, JA1225, JA1231, JA1235, JA1241 and JA1245

Raji (LDH-HiBiT) TCK Bioassays

All technical literature is available at: www.promega.com/protocols/
 Visit the website to verify that you are using the most current version of this Technical Manual.
 Email Promega Technical Services if you have questions on use of this system: techserv@promega.com

1. Description	3
2. Product Components and Storage Conditions	10
3. Before You Begin	13
4. PBMC ADCC Assay Protocol	13
4.A. Preparing Assay Reagents and Samples	16
4.B. Plate Layout Design	17
4.C. Thawing Human PBMC, ADCC-Qualified Effector Cells	18
4.D. Thawing Raji (LDH-HiBiT) Target Cells	18
4.E. Day 2: Preparing 3X Stock Solutions of Antibody	19
4.F. Day 2: Preparing 3X Digitonin for Maximum Release Control	19
4.G. Day 2: Preparing Raji Cells (LDH-HiBiT)	20
4.H. Day 2: Preparing Human PBMC, ADCC-Qualified Effector Cells	21
4.I. Day 2: Adding Target Cells, Effector Cells and Antibody	21
4.J. Day 2: Preparing and Adding Bio-Glo-NB™ TCK Reagent	21
4.K. Data Analysis	22
4.L. Using Knockout Controls	23
5. CD8+ T Cell TDCC Assay Protocol	23
5.A. Preparing Assay Reagents and Samples	25
5.B. Plate Layout Design	26
5.C. Thawing Raji Cells (LDH-HiBiT)	27
5.D. Day 2: Preparing 3X Stock Solutions of Biologic Drugs	28
5.E. Day 2: Preparing Raji Cells (LDH-HiBiT)	29
5.F. Day 2: Preparing Human CD8+ T Cells, TDCC Qualified Effector Cells (T Cells)	29
5.G. Day 2: Adding Raji Cells (LDH-HiBiT), Effector Cells and Biologic Drugs	29
5.H. Day 2: Preparing and Adding Bio-Glo-NB™ TCK Reagent	30
5.I. Data Analysis	31
5.J. Using Knockout Controls	31



- 6. CAR-T Assay Protocol 32
 - 6.A. Preparing Assay Reagents 34
 - 6.B. Plate Layout Design 35
 - 6.C. Thawing Raji Cells (LDH-HiBiT) 35
 - 6.D. Day 2: Preparing Raji Cells (LDH-HiBiT) 36
 - 6.E. Day 2: Preparing 2X Stock Solutions of CAR-T Cells 36
 - 6.F. Day 2: Adding Target Cells and CAR-T Cells to Assay Plates 37
 - 6.G. Day 2: Preparing and Adding Bio-Glo-NB™ TCK Reagent 37
 - 6.H. Data Analysis 38
 - 6.I. Using Knockout Controls 39
- 7. Troubleshooting 39
- 8. References 41
- 9. Appendix 42
 - 9.A. Composition of Buffers and Solutions 42
 - 9.B. Related Products 42

1. Description

Inducing tumor cell lysis by immune effector cells is a primary mechanism of action (MOA) for many immuno-oncology drugs. Monoclonal antibodies against tumor-associated antigens (TAA) can induce antibody-dependent cellular cytotoxicity (ADCC) by natural killer (NK) cells (1). Next-generation biologics, such as bispecific antibodies, can redirect T cells to lyse tumor cells (2). CAR-T cells and other cellular therapies are living drugs that directly lyse tumor cells as a component of their MOA (3). HiBiT Target Cell Killing (TCK) Bioassays provide a flexible, reproducible, label-free, luminescent system for measuring the potency of a variety of cellular and immuno-oncology drug products where target cell killing is a MOA.

Raji (LDH-HiBiT) TCK Bioassays^(a-7) measure the death of Raji cells engineered to express a fusion of HiBiT to the C terminus of the lactate dehydrogenase (LDHB) gene product (LDH-HiBiT). HiBiT is an 11 amino acid peptide that binds with high affinity to LgBiT to form NanoBiT[®] luciferase (4). When dead or dying cells release the HiBiT fusion protein into the extracellular medium, it binds to the cell-impermeable LgBiT, which is provided with the furimazine substrate in a nonlytic, homogeneous detection reagent (Figure 1). In coculture experiments, the luminescent signal will be proportional to target cell death alone, with no signal contribution arising from the death of effector cells.

Raji (LDH-HiBiT) TCK Bioassays are a suite of four individual bioassays comprised of the parental Raji line and three different knockout (KO) variants (see Section 2). Each bioassay contains a clonal cell line in the thaw-and-use format, which can be thawed, plated and used in assays without the need for continuous cell culture. All cell lines express the LDH-HiBiT reporter, making them suitable for analysis using the HiBiT TCK technology.


The Raji (LDH-HiBiT) TCK Bioassay has been prequalified in a T cell-dependent cellular cytotoxicity (TDCC) bioassay following International Council for Harmonisation of Technical Requirements for Pharmaceuticals for Human Use (ICH) guidelines and exhibit the precision, accuracy and linearity required for routine use in potency and stability studies (Table 1 and Figure 3). HiBiT TCK Bioassay workflows are simple, robust and compatible with both 96- and 384-well plate formats used for early biologic drug discovery and lot release settings (Figure 4).

Detection of released HiBiT fusion protein is sensitive, enabling the use of low numbers of HiBiT target cells (2,500–5,000 cells per well). The cells have a low rate of spontaneous HiBiT fusion protein release and the fusion protein is stable once released into the extracellular medium. These features enable use of the Raji (LDH-HiBiT) TCK Bioassays to measure the potency of biologic drugs and cell therapy products, as well as to monitor the specificity of drugs targeting CD19, CD20 or both receptors. For cell therapy, you can test a wide range of effector-to-target (E:T) ratios while conserving test material for experiments ranging from 4–72 hours (Figure 5).

When parental and knockout cells are cocultured with effector cells (e.g., NK cells, cytotoxic T cells, etc.) CD19-, CD20- or CD19/CD20-dependent killing can be monitored. For example, biologic drugs that target antigens expressed on the cell surface show dose-dependent signal increases (Figure 6). In contrast, when knockout cells are assayed with biologic drugs that target antigens no longer expressed on the cell surface, dose-dependent signal increases are absent (Figure 6).

This technical manual provides protocols for the suite of Raji (LDH-HiBiT) TCK Bioassays, including sample protocols for antibody-dependent cellular cytotoxicity (ADCC) bioassays using Human PBMC, ADCC-Qualified cells (Section 4); TDCC bioassays using Human T Cell (CD8+), TDCC-Qualified cells (Section 5); and CAR-T bioassays with user-provided CAR-T cells (Section 6). Optimization may be necessary for use with your biologic product of interest.

1. Description (continued)



193477A

Figure 1. Representation of the HiBiT TCK Bioassay. Target cells stably expressing a HiBiT fusion protein are cocultured with effector cells (CAR-T, T cells, NK cells, etc.). Effector cell-mediated killing of target cells releases the HiBiT fusion protein into the extracellular medium. Cell-impermeable LgBiT and furimazine substrate (Fz) are added as components of the Bio-Glo-NB™ TCK Reagent. HiBiT complementation with LgBiT generates NanoBiT® luciferase, a bright, luminescent enzyme. (Created with BioRender.com.)





Figure 2. Raji (LDH-HiBiT) TCK Bioassays reflect the MOA of biologic drugs and demonstrates specificity for cell surface antigens. Panel A. Raji Cells (LDH-HiBiT) were incubated with Human PBMC, ADCC-Qualified effector cells in the presence of serial antibody titrations, as indicated. After a 5-hour induction, Bio-Glo-NB™ TCK Reagent was added and luminescence quantified using the GloMax® Discover System. Raji Cells (LDH-HiBiT) express CD19 and CD20, but they do not express EGFR or HER2 on the cell surface. **Panel B.** Raji Cells (LDH-HiBiT) were incubated with Human T Cells (CD8+) in the presence of serial T cell engager (TCE) titrations, as indicated. After a 24-hour induction, Bio-Glo-NB™ TCK Reagent was added and luminescence quantified using the GloMax® Discover System. Raji Cells (LDH-HiBiT) express CD19 and BCMA, but they do not express HER2 or GPC3 on the cell surface. Data were fitted to a four-parameter logistic curve using GraphPad Prism® software. Data were generated using thaw-and-use cells.

1. Description (continued)

Table 1. Precision, Accuracy and Linearity of TDCC Bioassay with Raji Cells (LDH-HiBiT) as Targets.

Parameter	Results	
Accuracy	% Expected Relative Potency	% Recovery
	50	102.9
	75	104.5
	100	104.8
	125	108.5
	150	99.7
Repeatability (% CV)	100% (Reference)	
Intermediate Precision (% CV)	9.5	
Linearity (r ²)	0.93	
Linearity (y = mx + b)	y = 1.014x + 2.520	
<p>A 50–150% theoretical potency series of blinatumomab was analyzed in triplicate in three independent experiments performed on three days by two analysts using Raji Cells (LDH-HiBiT) and Human T Cells (CD8+), TDCC-Qualified effectors. Bio-Glo-NB™ TCK Reagent was added and luminescence quantified using the GloMax® Discover System. Data were analyzed and relative potencies calculated after parallelism determination using JMP® software. Data were generated using thaw-and-use cells.</p>		



189411AA

Figure 3. Raji (LDH-HiBiT) TCK Bioassays are stability-indicating. Panel A. Samples of rituximab were maintained at 4°C (control) or heat-treated at 65°C for the indicated times, then analyzed in an ADCC Bioassay with Raji Cells (LDH-HiBiT) and Human PBMC, ADCC-Qualified effector cells. **Panel B.** Samples of blinatumomab (CD19 × CD3) were maintained at 4°C (control) or heat-treated at 55°C for the indicated times, then analyzed in a TDCC Bioassay with Raji Cells (LDH-HiBiT) and Human T Cells (CD8+), TDCC-Qualified. For both panels, Bio-Glo-NB™ TCK Reagent was added and luminescence quantified using the GloMax® Discover System. Data were fitted to a four-parameter logistic curve using GraphPad Prism® software. Data were generated using thaw-and-use cells.

1. Description (continued)




Figure 4. Raji (LDH-HiBiT) TCK Bioassays are amenable to 384-well plate format. Panel A. Raji Cells (LDH-HiBiT) were harvested after overnight culture and 1×10^3 cells/5 μ l/well plated in a 384-well round-bottom white assay plate (e.g., Corning® Cat.# 4512, low volume). Next, 5 μ l of 3X serially-diluted rituximab was added, then 2.5×10^4 /5 μ l/well of Human PBMC, ADCC-Qualified effector cells from overnight culture were added. After a 5-hour incubation at 37°C, 5% CO₂, 5 μ l of 4X Bio-Glo-NB™ TCK Reagent was added to each well and luminescence quantified. **Panel B.** Raji Cells (LDH-HiBiT) were plated at 5×10^2 cells/5 μ l/well in a 384-well flat-bottom white assay plate (e.g., Corning® Cat.# 3570). Next, 5 μ l of 3X serially diluted blinatumomab was added, followed by 5×10^3 /5 μ l/well of Human T Cells (CD8+), TDCC-Qualified. After a 24-hour incubation at 37°C, 5% CO₂, 15 μ l of Bio-Glo-NB™ TCK Reagent was added per well and luminescence quantified. Luminescence was measured using the GloMax® Discover System. Data were fitted to four-parameter logistic curves using GraphPad Prism® software. Data were generated using thaw-and-use cells.




Figure 5. Raji (LDH-HiBiT) TCK Bioassays are amenable to use with CAR-T effector cells. **Panel A.** Raji Cells (LDH-HiBiT) were incubated with serially diluted human CD19 CAR-T cells at the indicated effector-to-target (E:T) ratios. After a 4-, 24-, 48- or 72-hour induction, Bio-Glo-NB™ TCK Reagent was added and luminescence quantified using the GloMax® Discover System. **Panel B.** The target cell killing data are represented as a percent of the maximum CAR-T-dependent lysis at each time point. Data were fitted to a four-parameter logistic curve using GraphPad Prism® software. Luminescence increased over time (**Panel A**) and EC₅₀ values decreased (data not shown) with longer incubation periods (**Panel B**). Data were generated using thaw-and-use cells.

1. Description (continued)




Figure 6. Raji (LDH-HiBiT) TCK Bioassays reflect biologic drug activity in wild-type and surface antigen knockout cells. Parental and knockout cells were incubated with Human T Cells (CD8+) in the presence of serial titrations of biologic drugs for 24 hours. **Panel A.** Raji cell lines expressing CD19 show a dose-dependent increase in percent specific lysis values for blinatumumab (CD19 × CD3), whereas CD19 KO cell lines do not. **Panel B.** Raji cell lines expressing CD20 show a dose-dependent increase in percent specific lysis values for mosunetuzumab (CD20 × CD3), whereas CD20 KO cell lines do not. **Panel C.** Raji cells express BCMA and show a dose-dependent increase in percent specific lysis values. Bio-Glo-NB™ TCK Reagent was added and luminescence measured using the GloMax® Discover System. Data were fitted to a four-parameter logistic curve using GraphPad Prism® software. Data were generated using thaw-and-use cells.

2. Product Components and Storage Conditions

PRODUCT	SIZE	CAT. #
Raji (LDH-HiBiT) TCK Bioassay	1 each	JA1311

Not for Medical Diagnostic Use. Each kit contains sufficient reagents for 120 assays using the inner 60 wells of two 96-well plates. Includes:

- 1 vial Raji Cells (LDH-HiBiT) (0.5ml per vial)
- 3 × 36ml RPMI 1640 Medium
- 2 × 4ml Fetal Bovine Serum
- 1 vial Bio-Glo-NB™ TCK Luciferase Assay Substrate
- 1 vial Bio-Glo-NB™ TCK LgBiT Protein
- 10ml Bio-Glo-NB™ TCK Luciferase Assay Buffer

PRODUCT	SIZE	CAT. #
Raji (LDH-HiBiT) TCK Bioassay, 5X	1 each	JA1315

Not for Medical Diagnostic Use. Each kit contains sufficient reagents for 600 assays using the inner 60 wells of two 96-well plates. Includes:

- 5 vials Raji Cells (LDH-HiBiT) (0.5ml per vial)
- 15 × 36ml RPMI 1640 Medium
- 10 × 4ml Fetal Bovine Serum
- 5 vials Bio-Glo-NB™ TCK Luciferase Assay Substrate
- 5 vials Bio-Glo-NB™ TCK LgBiT Protein
- 5 × 10ml Bio-Glo-NB™ TCK Luciferase Assay Buffer

PRODUCT	SIZE	CAT. #
Raji (LDH-HiBiT) CD19-KO TCK Bioassay	1 each	JA1221

Not for Medical Diagnostic Use. Each kit contains sufficient reagents for 120 assays using the inner 60 wells of two 96-well plates. Includes:

- 1 vial Raji Cells, CD19-KO (LDH-HiBiT) (0.5ml per vial)
- 3 × 36ml RPMI 1640 Medium
- 2 × 4ml Fetal Bovine Serum
- 1 vial Bio-Glo-NB™ TCK Luciferase Assay Substrate
- 1 vial Bio-Glo-NB™ TCK LgBiT Protein
- 10ml Bio-Glo-NB™ TCK Luciferase Assay Buffer

PRODUCT	SIZE	CAT. #
Raji (LDH-HiBiT) CD19-KO TCK Bioassay, 5X	1 each	JA1225

Not for Medical Diagnostic Use. Each kit contains sufficient reagents for 600 assays using the inner 60 wells of two 96-well plates. Includes:

- 5 vials Raji Cells, CD19-KO (LDH-HiBiT) (0.5ml per vial)
- 15 × 36ml RPMI 1640 Medium
- 10 × 4ml Fetal Bovine Serum
- 5 vials Bio-Glo-NB™ TCK Luciferase Assay Substrate
- 5 vials Bio-Glo-NB™ TCK LgBiT Protein
- 5 × 10ml Bio-Glo-NB™ TCK Luciferase Assay Buffer

PRODUCT	SIZE	CAT. #
Raji (LDH-HiBiT) CD20-KO TCK Bioassay	1 each	JA1231

Not for Medical Diagnostic Use. Each kit contains sufficient reagents for 120 assays using the inner 60 wells of two 96-well plates. Includes:

- 1 vial Raji Cells, CD20-KO (LDH-HiBiT) (0.5ml per vial)
- 3 × 36ml RPMI 1640 Medium
- 2 × 4ml Fetal Bovine Serum
- 1 vial Bio-Glo-NB™ TCK Luciferase Assay Substrate
- 1 vial Bio-Glo-NB™ TCK LgBiT Protein
- 10ml Bio-Glo-NB™ TCK Luciferase Assay Buffer



2. Product Components and Storage Conditions (continued)

PRODUCT	SIZE	CAT. #
Raji (LDH-HiBiT) CD20-KO TCK Bioassay, 5X	1 each	JA1235

Not for Medical Diagnostic Use. Each kit contains sufficient reagents for 600 assays using the inner 60 wells of two 96-well plates. Includes:

- 5 vials Raji Cells, CD20-KO (LDH-HiBiT) (0.5ml per vial)
- 15 × 36ml RPMI 1640 Medium
- 10 × 4ml Fetal Bovine Serum
- 5 vials Bio-Glo-NB™ TCK Luciferase Assay Substrate
- 5 vials Bio-Glo-NB™ TCK LgBiT Protein
- 5 × 10ml Bio-Glo-NB™ TCK Luciferase Assay Buffer

PRODUCT	SIZE	CAT. #
Raji (LDH-HiBiT) CD19/20-KO TCK Bioassay	1 each	JA1241

Not for Medical Diagnostic Use. Each kit contains sufficient reagents for 120 assays using the inner 60 wells of two 96-well plates. Includes:

- 1 vial Raji Cells, CD19/20-KO (LDH-HiBiT) (0.5ml per vial)
- 3 × 36ml RPMI 1640 Medium
- 2 × 4ml Fetal Bovine Serum
- 1 vial Bio-Glo-NB™ TCK Luciferase Assay Substrate
- 1 vial Bio-Glo-NB™ TCK LgBiT Protein
- 10ml Bio-Glo-NB™ TCK Luciferase Assay Buffer

PRODUCT	SIZE	CAT. #
Raji (LDH-HiBiT) CD19/20-KO TCK Bioassay, 5X	1 each	JA1245

Not for Medical Diagnostic Use. Each kit contains sufficient reagents for 600 assays using the inner 60 wells of two 96-well plates. Includes:


- 5 vials Raji Cells, CD19/20-KO (LDH-HiBiT) (0.5ml per vial)
- 15 × 36ml RPMI 1640 Medium
- 10 × 4ml Fetal Bovine Serum
- 5 vials Bio-Glo-NB™ TCK Luciferase Assay Substrate
- 5 vials Bio-Glo-NB™ TCK LgBiT Protein
- 5 × 10ml Bio-Glo-NB™ TCK Luciferase Assay Buffer

Storage Conditions: Upon arrival, immediately transfer the cell vials to below -140°C (freezer or liquid nitrogen vapor phase) for long-term storage. Do not store cell vials submerged in liquid nitrogen. Do not store cell vials at -80°C because this will decrease cell viability and cell performance. Store Bio-Glo-NB™ TCK Luciferase Assay Substrate, Bio-Glo-NB™ TCK LgBiT Protein, Bio-Glo-NB™ TCK Luciferase Assay Buffer and Fetal Bovine Serum (FBS) at -30°C to -10°C . Avoid multiple freeze-thaw cycles of the serum. The Bio-Glo-NB™ TCK Luciferase Assay Substrate and Bio-Glo-NB™-TCK LgBiT Protein remain as liquids and do not freeze. Store RPMI 1640 Medium at $+2^{\circ}\text{C}$ to $+10^{\circ}\text{C}$, protected from fluorescent light.

3. Before You Begin

Please read through the entire protocol to become familiar with the components and the assay procedures before beginning.

Remove the product label from the box containing vials of cells or note the catalog number, lot number and dispensed lot number from the label. This information can be used to download specific documents from www.promega.com, such as the Certificate of Analysis.

 **Note:** HiBiT TCK Bioassays use the Bio-Glo-NB™ TCK Luciferase Assay System (Cat.# JB1001, JB1002, JB1003) for detection. **Do not** use the Bio-Glo-NL™ Luciferase Assay System (Cat.# J3081, J3082, J3083) or the Bio-Glo™ Luciferase Assay System (Cat.# G7940, G7941).

Raji (LDH-HiBiT) TCK Bioassays are intended for use with effector cells capable of binding to and killing target cells alone or in combination with a biologic drug. Human PBMC, ADCC-Qualified, or Human T-Cell (CD8+), TDCC-Qualified, are available separately for use in ADCC and TDCC assays, respectively. For ADCC assays, Control Ab, Anti-CD20 (Cat.# GA1130) is available separately.

Raji (LDH-HiBiT) TCK Bioassays are provided in frozen, thaw-and-use format and are ready to be used without any additional cell culture or propagation. When thawed and diluted as instructed, the cells will be at the appropriate concentration for the assay. The cells are sensitive and care should be taken to follow cell thawing and plating procedures as described. Do not overmix or overwarm the cell reagents.

Raji (LDH-HiBiT) TCK Bioassays produce a bioluminescent signal and require a sensitive luminescence plate reader. Bioassay development and performance data included in this technical manual were generated using the GloMax® Discover System (see Section 9.C, Related Products). An integration time of 0.5 second/well was used for all readings.

Raji (LDH-HiBiT) TCK Bioassays are compatible with most other plate-reading luminometers, though relative light unit (RLU) readings will vary with the sensitivity and settings of each instrument. If using a reader with adjustable gain, we recommend a high-gain setting. The use of different instruments and gain adjustment will affect the magnitude of the raw data but should not affect the measured relative potency of test samples

4. PBMC ADCC Assay Protocol

This procedure illustrates the use of the Raji (LDH-HiBiT) TCK Bioassay to test two antibody samples against a reference sample in a single PBMC ADCC assay. Each test and reference antibody is run in triplicate, in an eight-point dilution series, in two 96-well assay plates using the inner 60 wells. Other experimental and plate layouts are possible but may require further optimization. See Section 4.L for recommendations on including knockout cell line controls.

4. PBMC ADCC Assay Protocol (continued)

Note: When preparing test and reference antibodies, choose an appropriate starting concentration and dilution scheme to achieve a complete dose-response curve with proper upper and lower asymptotes and sufficient points on the slope. For reference, we use 1 µg/ml as a starting concentration (1X) and threefold serial dilutions when testing Control Ab, Anti-CD20 (Cat.# GA1130).

Maximum release (MR) and spontaneous release (SR) controls should be included in your experiment. The MR control uses digitonin to permeabilize target cells, providing the maximal luminescent signal. The SR control is a measure of background cell death in presence of effector cells but in the absence of a monoclonal antibody.

Materials to Be Supplied By the User

(Composition of Buffers and Solutions is provided in Section 9.A.)

Reagents

- Human PBMC, ADCC-Qualified (e.g., www.promega.com/products/reporter-bioassays/primary-cell-bioassays/pbmc-adcc-bioassay/; please enquire)
- antibodies or other biologics capable of inducing ADCC
- recombinant human IL-2 (e.g., Sigma Cat.# I2644)
- β-mercaptoethanol (e.g., GIBCO® Cat.# 21985)
- sodium pyruvate (e.g., GIBCO® Cat.# 11360)
- digitonin (e.g., Cat.# G9441)
- Trypan blue solution (e.g., Sigma Cat.# T8154)
- **optional:** Control Ab, Anti-CD20 (Cat.# GA1130)

Supplies and Equipment

- solid-white, U-bottom 96-well assay plates (e.g., Corning® Cat.# 3355) or 384-well assay plates (e.g., Corning® Cat.# 4512, 4513) for ADCC assay applications
- sterile clear V-bottom 96-well plate with lid (e.g., Costar® Cat.# 3896) for preparing antibody dilutions
- pipettes (single-channel and 12-channel; for best results use both manual and electronic pipettes as needed)
- sterile 15ml and 50ml conical tubes
- sterile reagent reservoirs (e.g., Corning®/Costar® Cat.# 4870)
- T75 cell culture flasks (e.g., Corning® Cat.# 430641U)
- T25 cell culture flasks (e.g., Corning® Cat.# 430639)
- 37°C, 5% CO₂ incubator
- 37°C water bath
- sensitive plate reader with luminescence measuring capability or luminometer (e.g., GloMax® Discover System or equivalent)

Note: Corning® Cat.# 3355 plates are supplied without a lid and are marked as nonsterile. Use a lid from another Corning® plate type (e.g., Cat.# 3917). Although marked as nonsterile, we had no issues with sterility when using these plates while developing this product.




Figure 7. Schematic protocol for the Raji (LDH-HiBiT) TCK Bioassay in a PBMC ADCC assay application.

4.A. Preparing Assay Reagents and Samples

Complete formulas for the following reagents are provided in Section 9.A.

Heat-Inactivated (HI) FBS: Thaw the provided FBS overnight at +2°C to +10°C or in a 37°C water bath on the day of use. Heat inactivate the FBS for 30 minutes at 56°C. Mix well and adjust to 37°C before use.


Raji Cells (LDH-HiBiT) Recovery Medium: On the day before the assay, prepare 15ml of Raji Cells (LDH-HiBiT) recovery medium (see Section 9.A). Thaw both vials of FBS overnight at +2°C to +10°C or in a 37°C water bath on the day of use. Heat inactivate the FBS as noted above. Mix well and adjust to 37°C before use. To prepare Raji Cells (LDH-HiBiT) recovery medium, add 1.5ml of FBS to 13.5ml of 37°C RPMI 1640 Medium. For reference, 15ml of Raji Cells (LDH-HiBiT) recovery medium is sufficient to thaw and plate one vial of Raji Cells (LDH-HiBiT). If multiple vials will be thawed, then scale the amount of Raji Cells (LDH-HiBiT) recovery medium appropriately. Store the remaining FBS at +2°C to +10°C for use in preparing the assay buffer on the day of the assay.

PBMC Culture Medium: On the day before the assay, prepare an appropriate amount of PBMC culture medium: 90% RPMI 1640 with 10% HI-FBS and IL-2, β -mercaptoethanol, sodium pyruvate.

PBMC ADCC Assay Buffer: On the day of the assay, prepare an appropriate amount of assay buffer: 95% RPMI 1640 with 5% HI-FBS with IL-2 and β -mercaptoethanol. Mix well and warm to 37°C before use.

Note: The recommended assay buffer contains 5% HI-FBS. This concentration of FBS works well for the Control Ab, Anti-CD20 (Cat.# GA1130) that we tested. If you experience assay performance issues when using this assay buffer, we recommend testing different serum concentrations in the range of 0.5–10%.

Bio-Glo-NB™ TCK Reagent: For reference, 10ml of Bio-Glo-NB™ TCK Reagent is sufficient to assay 120 wells in a 96-well assay format. Store the Bio-Glo-NB™ TCK Luciferase Assay Substrate and LgBiT Protein at –30°C to –10°C. Thaw the Bio-Glo-NB™ TCK Luciferase Assay Buffer at room temperature (do not exceed 25°C) during the 5-hour assay induction period. We recommend preparing the reconstituted Bio-Glo-NB™ TCK Reagent immediately before use.

 **Note:** HiBiT TCK Bioassays are compatible only with Bio-Glo-NB™ TCK Reagent. **Do not** use Bio-Glo™ Reagent or Bio-Glo-NL™ Reagent.

Test and Reference Samples: Using assay buffer as the diluent, prepare stock starting dilutions (dilu1, 3X final concentration) of two test antibodies (300 μ l each) and one reference antibody (600 μ l) in tubes. Store the tubes containing antibody starting dilutions appropriately before making antibody serial dilutions.

Note: If you are using Control Ab, Anti-CD20 (Cat.# GA1130), as a reference antibody in your assay, prepare 300 μ l of a starting dilution at 3 μ g/ml (dilu1, 3X final concentration) in PBMC ADCC assay buffer.


4.B. Plate Layout Design

For the protocol described here, use the plate layout illustrated in Figure 8 as a guide. The protocol describes serial replicate dilutions (n = 3) of test and reference antibodies to generate two eight-point dose-response curves for each plate.

Recommended Plate Layout Design													
	1	2	3	4	5	6	7	8	9	10	11	12	
A	B	B	B	B	B	B	B	B	B	B	B	B	Assay Buffer (B)
B	B	SR	dilu8	dilu7	dilu6	dilu5	dilu4	dilu3	dilu2	dilu1	MR	B	Reference Ab
C	B	SR	dilu8	dilu7	dilu6	dilu5	dilu4	dilu3	dilu2	dilu1	MR	B	Test Ab
D	B	SR	dilu8	dilu7	dilu6	dilu5	dilu4	dilu3	dilu2	dilu1	MR	B	Reference Ab
E	B	SR	dilu8	dilu7	dilu6	dilu5	dilu4	dilu3	dilu2	dilu1	MR	B	Test Ab
F	B	SR	dilu8	dilu7	dilu6	dilu5	dilu4	dilu3	dilu2	dilu1	MR	B	Reference Ab
G	B	SR	dilu8	dilu7	dilu6	dilu5	dilu4	dilu3	dilu2	dilu1	MR	B	Test Ab
H	B	B	B	B	B	B	B	B	B	B	B	B	Assay Buffer (B)

Figure 8. Example plate layout showing nonclustered sample locations of test antibody and reference antibody dilution series. Maximum release (MR) control wells contain Raji cells (LDH-HiBiT), PBMC effector cells and digitonin. Spontaneous release (SR) control wells contain Raji Cells (LDH-HiBiT) and PBMC effector cells. Test or reference antibodies are not added to MR and SR control wells. Wells containing assay buffer alone are denoted by "B".


4.C. Thawing Human PBMC, ADCC-Qualified Effector Cells

 Follow institutional guidelines for handling, including use of personal protective equipment (PPE) and waste disposal for biohazardous material.

Note: Perform the following steps in a sterile cell culture hood.

1. Remove one vial ($\sim 2 \times 10^7$) of PBMC from storage at -140°C and thaw in a 37°C water bath with gentle agitation (do not invert cell vial) until just thawed (typically 2–3 minutes).
2. Transfer all of the cells (approximately 1ml) to a 50ml conical tube containing 24ml of prewarmed PBMC ADCC culture medium.
3. Gently resuspend the PBMC using a 25ml pipette.
4. Transfer the cell suspension to a T75 tissue culture flask and place the flask horizontally in a humidified 37°C , 5% CO_2 incubator.
5. Incubate for approximately 16–24 hours before assay setup.

4.D. Thawing Raji (LDH-HiBiT) Target Cells

 **Note:** Perform the following steps in a sterile cell culture hood.

Thaw-and-use Raji Cells (LDH-HiBiT) are in a sensitive state immediately post-thaw and care should be taken to follow the thawing and plating procedures exactly as described. Handle the cells gently and do not overmix the cells or overwarm the reagents. Do not vortex the cells or centrifuge at speeds $>150 \times g$. It is important to keep pipette tips away from tube surfaces during delivery or mixing steps to minimize mechanical stress on the cells. We recommend that you thaw and dilute a maximum of two vials of thaw-and-use cells at any one time.

1. Add 9.5ml of Raji Cells (LDH-HiBiT) recovery medium to a 15ml conical tube and warm to 37°C .
2. Remove one cryovial of Raji Cells (LDH-HiBiT) from storage at -140°C and thaw in a 37°C water bath with gentle agitation (do not invert) until just thawed (typically 2–3 minutes).
3. Using a 2ml serological pipette or a 1ml wide-bore pipette tip, transfer 1ml of prewarmed Raji Cells (LDH-HiBiT) recovery medium, one drop at a time, to the cryovial. Using the same pipette, gently mix the contents of the cryovial by slowly pipetting approximately 1ml, once. Gently transfer the contents of the cryovial to the 15ml conical tube containing prewarmed Raji Cells (LDH-HiBiT) recovery medium.
4. Pellet the cells at $150 \times g$ for 5 minutes.
5. Aspirate the supernatant, leaving 100–300 μl of Raji Cells (LDH-HiBiT) recovery medium above the cell pellet.
6. Resuspend the cell pellet by gently flicking the tube bottom 2–3 times. Add 5ml of prewarmed Raji Cells (LDH-HiBiT) recovery medium to the side of the tube. Do not pipet Raji Cells (LDH-HiBiT) recovery medium directly onto the cell pellet. Pipet once to mix.
7. Transfer the cell suspension to a T25 cell culture flask and place the flask horizontally in a humidified 37°C , 5% CO_2 incubator.
8. Incubate for approximately 16–24 hours before assay setup.

4.E. Day 2: Preparing 3X Stock Solutions of Antibody

The instructions described here are for the preparation of 3X stocks of test and reference antibody. A threefold dilution series is made for each antibody, where a single 150µl stock of each dilution provides sufficient volume for analysis in triplicate. To prepare threefold serial dilutions, you will need 300µl of each test antibody and 600µl of reference antibody at 3X the highest antibody concentration in your dose-response curve. For other dilution schemes, adjust the volumes accordingly. For instance, prepare three independent stocks of each dilution for analysis in triplicate.

1. On the day of the assay, prepare an appropriate amount of PBMC ADCC assay buffer as described in Section 9.A.
2. To a sterile clear V-bottom 96-well plate, add 300µl of reference antibody starting dilution (dilu1, 3X final concentration) to wells A10 and B10 (Figure 9).
3. Add 300µl of test antibodies 1 and 2 starting dilution (dilu1, 3X final concentration) to wells E10 and G10, respectively (Figure 9).
4. Add 200µl of PBMC ADCC assay buffer to the other wells in these four rows, from column 9 to column 2.
5. Transfer 100µl of the antibody starting dilutions from column 10 into column 9. Mix well by pipetting. Avoid creating bubbles.
6. Repeat equivalent threefold serial dilutions across the columns from right to left through column 3. Do not dilute into column 2.
Note: Wells A2, B2, E2 and G2 contain 200µl of PBMC ADCC assay buffer without antibody as a negative control for delivery to SR control wells.
7. Proceed immediately to Section 4.F.

4.F. Day 2: Preparing 3X Digitonin for Maximum Release Control

1. Dilute digitonin stock in PBMC ADCC assay buffer to 300µg/ml (3X). Prepare 100µl for each test antibody and 200µl for the reference antibody.
2. Transfer 100µl of the digitonin solution to wells A11, B11, E11 and G11 of the antibody dilution plate.
3. Cover the antibody dilution plate with a lid and incubate at 37°C in a humidified 5% CO₂ incubator until Section 4.I.

4.F. Day 2: Preparing 3X Digitonin for Maximum Release Control (continued)

Recommended Plate Layout for Antibody Dilutions Prepared from a Single Antibody Stock.													
	1	2	3	4	5	6	7	8	9	10	11	12	
A		no Ab	dilu8	dilu7	dilu6	dilu5	dilu4	dilu3	dilu2	dilu1	digitonin		Reference Ab
B		no Ab	dilu8	dilu7	dilu6	dilu5	dilu4	dilu3	dilu2	dilu1	digitonin		Reference Ab
C													
D													
E		no Ab	dilu8	dilu7	dilu6	dilu5	dilu4	dilu3	dilu2	dilu1	digitonin		Test Ab 1
F													
G		no Ab	dilu8	dilu7	dilu6	dilu5	dilu4	dilu3	dilu2	dilu1	digitonin		Test Ab 2
H													

Figure 9. Antibody dilution plate. Test and reference antibodies are serially diluted prior to addition to assay plates. Digitonin-containing samples in column 11 are added to MR control wells. No antibody samples in column 2 are added to SR control wells.

4.G. Day 2: Preparing Raji Cells (LDH-HiBiT)



Note: Perform the following steps in a sterile cell culture hood.

1. Gently mix and harvest the Raji Cells (LDH-HiBiT) that were incubated overnight.
2. Transfer into a 15ml centrifuge tube.
3. Centrifuge cells at 150 × g for 5 minutes.
4. Wash cells by removing supernatant, and gently resuspending in 10ml of prewarmed assay buffer.
5. Centrifuge cells at 150 × g for 5 minutes, then remove and discard supernatant.
6. Gently resuspend the cell pellet in 10ml of prewarmed assay buffer.
7. Incubate the Raji Cells (LDH-HiBiT) at ambient temperature until Section 4.I.

4.H. Day 2: Preparing Human PBMC, ADCC-Qualified Effector Cells

1. Gently resuspend and harvest PBMC using a 25ml pipette.
2. Transfer all of the PBMC to a 50ml conical centrifuge tube.
3. Pellet the cells at $300 \times g$ for 5 minutes at ambient temperature.
4. Gently but thoroughly resuspend the pellet in 4ml of PBMC ADCC assay buffer to a cell density of 5×10^6 cells/ml.
5. Keep the PBMC at ambient temperature until Section 4.I.

4.I. Day 2: Adding Target Cells, Effector Cells and Antibody

 **Note:** Use only white assay plates with U-shaped well bottoms.


1. Gently resuspend and transfer the Raji Cells (LDH-HiBiT) suspension to a sterile reagent reservoir.
2. Using a multichannel pipette, immediately dispense 25 μ l of the cell suspension to wells B2–G11 of a 96-well white U-bottom assay plate. The final cell number in each well should be 5×10^3 cells/well.
3. Using a multichannel pipette, add 25 μ l of the appropriate antibody dilution to wells B3–G10 of the assay plate according to Figure 8. Add 25 μ l of PBMC ADCC assay buffer to SR control wells (B2–G2). Add 25 μ l of PBMC ADCC assay buffer plus 300 μ g/ml digitonin to MR control wells (B11–G11).
4. Gently resuspend and transfer the PBMC effector cell suspension to a sterile reagent reservoir.
5. Using a multichannel pipette, immediately dispense 25 μ l of the cell suspension to each of wells B2–G11 of the assay plates. The final cell number in each well should be 1.25×10^5 cells/well.
6. Add 75 μ l of PBMC ADCC assay buffer to each of the outside wells of the assay plates.
7. Cover the assay plates with lids and incubate in a 37°C, 5% CO₂ incubator for 5 hours.

Note: The 5-hour incubation was optimized using Control Ab, Anti-CD20. We recommend optimizing assay time (3–24 hours) with your antibody or other biologic samples.

8. During assay incubation, warm Bio-Glo-NB™ TCK Luciferase Assay Buffer to ambient temperature.

4.J. Day 2: Preparing and Adding Bio-Glo-NB™ TCK Reagent

We recommend preparing the Bio-Glo-NB™ TCK Reagent immediately before use. Ensure that Bio-Glo-NB™ TCK Luciferase Assay Buffer is equilibrated to room temperature (do not exceed 25°C) before reconstituting the reagent. Do not store the reconstituted reagent. Once reconstituted, the reagent will lose about 15% activity over 8 hours and about 60% activity over 24 hours at room temperature.

 **Note:** HiBiT TCK Bioassays are compatible only with Bio-Glo-NB™ TCK Reagent. **Do not** use Bio-Glo™ Reagent or Bio-Glo-NL™ Reagent.

1. Remove the Bio-Glo-NB™ TCK Luciferase Assay Substrate from –30°C to –10°C storage and mix by pipetting. Briefly centrifuge the tubes if the substrate has collected in the cap or on the sides of the tubes.

4.J. Day 2: Preparing and Adding Bio-Glo-NB™ TCK Reagent (continued)

2. Remove the Bio-Glo-NB™ TCK LgBiT Protein from –30°C to –10°C storage and mix by pipetting. Briefly centrifuge the tubes if the solution has collected in the cap or on the sides of the tubes.
3. Transfer the desired amount of ambient temperature Bio-Glo-NB™ TCK Luciferase Assay Buffer to a 15ml or 50ml centrifuge tube.
4. Add Bio-Glo-NB™ TCK LgBiT Protein (1:100 dilution) and Bio-Glo-NB™ TCK Luciferase Assay Substrate (1:50 dilution) to the Bio-Glo-NB™ TCK Luciferase Assay Buffer. For example, if the experiment requires 10ml of reagent, add 100µl of Bio-Glo-NB™ TCK LgBiT Protein and 200µl of Bio-Glo-NB™ TCK Luciferase Assay Substrate to 10ml of Bio-Glo-NB™ TCK Luciferase Assay Buffer. Ten milliliters (10ml) of Bio-Glo-NB™ TCK Reagent is sufficient for 120 wells (two assay plates, using the inner 60 wells of each plate).
5. Remove assay plates from the incubator after the incubation period and equilibrate to room temperature for 10–15 minutes.
6. Using a multichannel pipette, add 75µl of Bio-Glo-NB™ TCK Reagent to all assay wells, taking care to not create bubbles. This includes MR and SR control wells.
7. Wait 10 minutes, then measure the luminescence in a GloMax® Discover System or a plate reader with glow-type luminescence reading capabilities. The luminescence intensity will decay gradually, with a signal half-life of 1–2 hours at room temperature.

Note: Varying the Bio-Glo-NB™ TCK incubation time will affect the raw relative light unit (RLU) values but should not significantly change the EC₅₀ value and maximum fold induction.

4.K. Data Analysis

The HiBiT TCK Bioassay is a nonlytic bioassay that measures HiBiT released from dead or dying target cells using Bio-Glo-NB™ TCK Reagent. Luminescence signal is expressed as RLU.

1. Calculate fold induction:

$$\text{Fold Induction} = \frac{\text{RLU}_{\text{induced}}}{\text{RLU}_{\text{SR control}}}$$

2. Calculate percent specific lysis:

$$\text{Percent Specific Lysis} = \frac{\text{RLU}_{\text{induced}} - \text{RLU}_{\text{SR control}}}{\text{RLU}_{\text{MR control}} - \text{RLU}_{\text{SR control}}} \times 100$$

MR (maximum release) control wells contain Raji Cells (LDH-HiBiT) plus PBMC effector cells (no test or reference antibody added). These wells are treated with digitonin to measure the total amount of HiBiT protein in each well. SR (spontaneous release) control wells contain Raji Cells (LDH-HiBiT) plus PBMC effector cells (no test or reference antibody added). These wells measure background killing in the absence of a test or reference antibody.

3. Graph data as RLU versus Log_{10} [antibody], fold induction versus Log_{10} [antibody] and percent specific lysis versus Log_{10} [antibody]. Fit curves and determine the EC_{50} value of antibody response using appropriate curve fitting software (such as GraphPad Prism[®] software).

4.L. Using Knockout Controls

Knockout cells should be plated at the same density as the parental cell line and handled identically during bioassay setup. You should continue to include test and reference antibodies on the same plate, using separate plates for parental and knockout cell lines. An antibody that targets an antigen remaining on the surface of knockout cells (e.g., CD22 or CD38) can be included as a positive control for both parental and knockout cell lines. Always calculate percent specific lysis when comparing between parental and knockout controls, as clones may vary in raw RLU.

5. CD8+ T Cell TDCC Assay Protocol

This procedure illustrates the use of the Raji (LDH-HiBiT) TCK Bioassay to test two samples against a reference sample in a single TDCC assay. Each test and reference sample is run in triplicate, in an eight-point dilution series, in a single 96-well assay plate using the inner 60 wells (Figure 10). Other experimental and plate layouts are possible but may require further optimization. See Section 5.J for recommendations on including knockout cell line controls.

Note: When preparing test and reference samples, choose an appropriate starting concentration and dilution scheme to achieve a complete dose-response curve with proper upper and lower asymptotes and sufficient points on the slope. For reference, we use 15ng/ml as a starting concentration (1X) and threefold serial dilutions when testing blinatumomab.

Maximum release (MR) and spontaneous release (SR) controls should be included in your experiment. The MR control uses digitonin to permeabilize target cells, providing the maximal luminescent signal. The SR control is a measure of background cell death in the absence of a biologic drug.

Materials to Be Supplied By the User

(Composition of Buffers and Solutions is provided in Section 9.A.)

Reagents

- Human T Cells (CD8+), TDCC-Qualified (please enquire: EarlyAccess@promega.com)
- bispecific antibodies or related biologic drugs capable of inducing TDCC
- digitonin (e.g., Cat.# G9441)

Supplies and Equipment

- solid-white, flat-bottom 96-well assay plates (e.g., Corning[®] Cat.# 3917) or 384-well assay plates (e.g., Corning[®] Cat.# 4512, 4513)
- sterile clear V-bottom 96-well plate with lid (e.g., Costar[®] Cat.# 3896) for preparing antibody dilutions
- pipettes (single-channel and 12-channel; for best results use both manual and electronic pipettes as needed)
- sterile 15ml and 50ml conical tubes
- sterile reagent reservoirs (e.g., Corning[®]/Costar[®] Cat.# 4870)
- T25 cell culture flask (e.g., Corning[®] Cat.# 430639)
- 37°C, 5% CO₂ incubator
- 37°C water bath
- sensitive plate reader with luminescence measuring capability or luminometer (e.g., GloMax[®] Discover System [Cat.# GM3000] or equivalent)

5. CD8+ T Cell TDC Assay Protocol (continued)




Figure 10. Schematic protocol for the Raji (LDH-HiBiT) TCK Bioassay in a TDC bioassay application.

5.A. Preparing Assay Reagents and Samples


Heat-Inactivated (HI) FBS: Thaw the provided FBS overnight at +2°C to +10°C or in a 37°C water bath on the day of use. Heat inactivate the FBS for 30 minutes at 56°C. Mix well and adjust to 37°C before use.

Raji Cells (LDH-HiBiT) Recovery Medium: On the day before the assay, prepare 15ml of Raji Cells (LDH-HiBiT) recovery medium (see Section 9.A). Thaw both vials of FBS overnight at +2°C to +10°C or in a 37°C water bath and heat inactivate as noted above. To prepare Raji Cells (LDH-HiBiT) recovery medium, add 1.5ml of FBS to 13.5ml of 37°C RPMI 1640 Medium. For reference, 15ml of Raji Cells (LDH-HiBiT) recovery medium is sufficient to thaw and plate one vial of Raji cells (LDH-HiBiT). If multiple vials will be thawed, then scale the amount of Raji Cells (LDH-HiBiT) recovery medium appropriately. Store the remaining HI-FBS at +2°C to +10°C for use in preparing the assay buffer on the day of the assay.

CD8+ T Cell TDCC Assay Buffer: On the day of the assay, prepare an appropriate amount of assay buffer (90% RPMI 1640, 10% HI-FBS). Mix well and warm to 37°C before use.

Note: The recommended assay buffer contains 10% HI-FBS. This concentration of FBS works well for use with blinatumomab. If you experience assay performance issues when using this assay buffer, we recommend testing serum concentrations in the range of 0.5–10%.

Bio-Glo-NB™ TCK Reagent: For reference, 10ml of Bio-Glo-NB™ TCK Reagent is sufficient to assay 120 wells in a 96-well assay format. The Bio-Glo-NB™ TCK Luciferase Assay Substrate and the Bio-Glo-NB™ TCK LgBiT Protein should always be stored at –30°C to –10°C. Thaw the Bio-Glo-NB™ TCK Luciferase Assay Buffer at room temperature (do not exceed 25°C) during the overnight assay induction period. We recommend preparing the reconstituted Bio-Glo-NB™ TCK Reagent immediately before use.

 **Note:** HIBiT TCK Bioassays are compatible only with Bio-Glo-NB™ TCK Reagent. **Do not** use Bio-Glo™ Reagent or Bio-Glo-NL™ Reagent.

Test and Reference Samples: Using assay buffer as the diluent, prepare stock starting dilutions (dilu1, 3X final concentration) of two test biologic drugs (300µl each) and one reference biologic (600µl) in tubes. Store the tubes containing antibody starting dilutions appropriately before making antibody serial dilutions.

5.B. Plate Layout Design

For the protocol described here, use the plate layout illustrated in Figure 11 as a guide. The protocol describes serial replicate dilutions (n = 3) of test and reference biologics to generate two eight-point dose-response curves for each plate.

Recommended Plate Layout Design													
	1	2	3	4	5	6	7	8	9	10	11	12	
A	B	B	B	B	B	B	B	B	B	B	B	B	Assay Buffer (B)
B	B	SR	dilu8	dilu7	dilu6	dilu5	dilu4	dilu3	dilu2	dilu1	MR	B	Reference Biologic
C	B	SR	dilu8	dilu7	dilu6	dilu5	dilu4	dilu3	dilu2	dilu1	MR	B	Test Biologic
D	B	SR	dilu8	dilu7	dilu6	dilu5	dilu4	dilu3	dilu2	dilu1	MR	B	Reference Biologic
E	B	SR	dilu8	dilu7	dilu6	dilu5	dilu4	dilu3	dilu2	dilu1	MR	B	Test Biologic
F	B	SR	dilu8	dilu7	dilu6	dilu5	dilu4	dilu3	dilu2	dilu1	MR	B	Reference Biologic
G	B	SR	dilu8	dilu7	dilu6	dilu5	dilu4	dilu3	dilu2	dilu1	MR	B	Test Biologic
H	B	B	B	B	B	B	B	B	B	B	B	B	Assay Buffer (B)

Figure 11. Example plate layout showing nonclustered sample locations of test and reference dilution series of biologic drugs. Maximum release (MR) and spontaneous release (SR) control wells contain Raji Cells (LDH-HiBiT) and CD8+ T cells but lack test or reference biologic drugs. MR control wells receive Bio-Glo NB™ TCK Reagent containing 200µg/ml digitonin. Wells containing assay buffer alone are denoted by “B”.

5.C. Thawing Raji Cells (LDH-HiBiT)

! **Note:** Perform the following steps in a sterile cell culture hood.

Thaw-and-use Raji Cells (LDH-HiBiT) are in a sensitive state immediately post-thaw and care should be taken to follow the thawing and plating procedures exactly as described. Handle the cells gently and do not overmix the cells or overwarm the reagents. Do not vortex the cells or centrifuge at speeds $>150 \times g$. It is important to keep pipette tips away from tube surfaces during delivery or mixing steps to minimize mechanical stress on the cells. We recommend that you thaw and dilute a maximum of two vials of thaw-and-use cells at any one time.

1. Add 9.5ml of Raji Cells (LDH-HiBiT) recovery medium to a 15ml conical tube and warm to 37°C.
2. Remove one cryovial of Raji Cells (LDH-HiBiT) from storage at -140°C and thaw in a 37°C water bath with gentle agitation (do not invert) until just thawed (typically 2–3 minutes).
3. Using a 2ml serological pipette or a 1ml wide-bore pipette tip, transfer 1ml of prewarmed Raji Cells (LDH-HiBiT) recovery medium, one drop at a time, to the cryovial. Using the same pipette, gently mix the contents of the cryovial by slowly pipetting approximately 1ml, once. Gently transfer the contents of the cryovial to the 15ml conical tube containing prewarmed Raji Cells (LDH-HiBiT) recovery medium.
4. Pellet the cells at $150 \times g$ for 5 minutes.
5. Aspirate the supernatant, leaving 100–300 μl of Raji Cells (LDH-HiBiT) recovery medium above the cell pellet.
6. Resuspend the cell pellet by gently flicking the tube bottom 2–3 times. Add 5ml of prewarmed Raji Cells (LDH-HiBiT) recovery medium to the side of the tube. Do not pipet Raji Cells (LDH-HiBiT) recovery medium directly onto the cell pellet. Pipet once to mix.
7. Transfer the cell suspension to a T25 cell culture flask and place the flask horizontally in a humidified 37°C, 5% CO_2 incubator.
8. Incubate for approximately 16–24 hours before assay setup.

5.D. Day 2: Preparing 3X Stock Solutions of Biologic Drugs

The instructions described here are for preparation of 3X stocks of test and reference biologic drugs. A threefold dilution series is made for each biologic drug, where a single 150µl stock of each dilution provides sufficient volume for analysis in triplicate. To prepare threefold serial dilutions, you will need 300µl of test biologic drugs and 600µl of a reference biologic drug at 3X the highest concentration used in your dose-response curve. For other dilution schemes, adjust the volumes accordingly. For instance, prepare three independent stocks of each dilution for analysis in triplicate.

Recommended Plate Layout for Biologic Dilutions Prepared from a Single Antibody Stock.													
	1	2	3	4	5	6	7	8	9	10	11	12	
A		No Biologic	dilu8	dilu7	dilu6	dilu5	dilu4	dilu3	dilu2	dilu1	No Biologic		Reference Biologic
B		No Biologic	dilu8	dilu7	dilu6	dilu5	dilu4	dilu3	dilu2	dilu1	No Biologic		Reference Biologic
C													
D													
E		No Biologic	dilu8	dilu7	dilu6	dilu5	dilu4	dilu3	dilu2	dilu1	No Biologic		Test Biologic 1
F													
G		No Biologic	dilu8	dilu7	dilu6	dilu5	dilu4	dilu3	dilu2	dilu1	No Biologic		Test Biologic 2
H													


Figure 12. Biologic drug dilution plate. Test and reference biologic drugs are serially diluted prior to addition to assay plates. No biologic drug samples in columns 2 and 11 are added to SR and MR control wells, respectively. Do not add digitonin to the wells in column 11. Digitonin is added to MR control wells at the time of luminescence measurement (Section 5.H).

Note: Perform the following steps in a sterile cell culture hood.

1. On the day of the assay, prepare an appropriate amount of TDCC assay buffer as described in Section 9.A.
2. To a sterile clear V-bottom 96-well plate, add 300µl of reference biologic starting dilution (dilu1, 3X final concentration) to wells A10 and B10 (Figure 12).
3. Add 300µl of test biologics 1 and 2 starting dilution (dilu1, 3X final concentration) to wells E10 and G10, respectively (Figure 12).
4. Add 200µl of TDCC assay buffer to dilu2 through dilu8 in rows A, B, E and G, from column 9 to column 3. Add 200µl of TDCC assay buffer to no biologic drug controls wells (columns 2 and 11).

5. Transfer 100µl of the biologic drug starting dilutions from column 10 into column 9. Mix well by pipetting. Avoid creating bubbles.
6. Repeat equivalent threefold serial dilutions across the columns from right to left through column 3. Do not dilute into column 2.
Note: Wells A2, B2, E2 and G2 contain 200µl of TDCC assay buffer without biologic drug for delivery to SR control wells. Wells A11, B11, E11 and G11 contain 200µl of TDCC assay buffer without biologic drug or digitonin for delivery to MR control wells. Digitonin is added to MR control wells at the time of luminescence measurement (Section 5.H).
7. Cover the biologic dilution plate with a lid and incubate at 37°C in a humidified 5% CO₂ incubator while preparing the Raji Cells (LDH-HiBiT) and CD8+ T cells.

5.E. Day 2: Preparing Raji Cells (LDH-HiBiT)

 **Note:** Perform the following steps in a sterile cell culture hood.

1. Gently mix and harvest the Raji Cells (LDH-HiBiT) that were cultured overnight (Section 5.C).
2. Transfer into a 50ml centrifuge tube.
3. Centrifuge cells at 150 × g for 5 minutes.
4. Wash cells by removing supernatant, and gently resuspending in 10ml of prewarmed assay buffer.
5. Centrifuge cells at 150 × g for 5 minutes, then remove and discard supernatant.
6. Gently resuspend the cell pellet in 20ml of prewarmed assay buffer.
7. Incubate the Raji Cells (LDH-HiBiT) at ambient temperature until Section 5.G.

5.F. Day 2: Preparing Human CD8+ T Cells, TDCC Qualified Effector Cells (T Cells)

1. Remove one vial of T cells from storage at -140°C and thaw in a 37°C water bath with gentle agitation (do not invert cell vial) until just thawed (typically 2–3 minutes).
2. Transfer all of the cells into 3.5ml of prewarmed TDCC assay buffer for a concentration of 1 × 10⁶ cells/ml.
3. Incubate the T cells at ambient temperature until Section 5.G.

5.G. Day 2: Adding Raji Cells (LDH-HiBiT), Effector Cells and Biologic Drugs

1. Gently resuspend and transfer the Raji Cells (LDH-HiBiT) suspension to a sterile reagent reservoir.
2. Using a multichannel pipette, immediately dispense 25µl of the cell suspension to each of wells B2–G11 of a 96-well, white, flat-bottom assay plate. The final cell number in each well should be 2.5 × 10³ cells/well.
3. Using a multichannel pipette, add 25µl of the appropriate biologic drug dilution (Figure 12) to wells B3–G10 of the assay plate according to the plate layout in Figure 11. Add 25µl of TDCC assay buffer to MR and SR control wells (B2–G2 and B11–G11, respectively).

Note: Digitonin will be added to MR control wells just prior to measuring luminescence (Section 5.H).

4. Gently resuspend and transfer the T cell suspension to a sterile reagent reservoir.

5.G. Day 2: Adding Raji Cells (LDH-HiBiT), Effector Cells and Biologic Drugs (continued)


5. Using a multichannel pipette, immediately dispense 25µl of the T cell suspension to each of wells B2–G11 of the assay plates. The final cell number in each well should be 2.5×10^4 cells/well.
6. Add 75µl of TDCC assay buffer to each of the outside wells of the assay plates.
7. Cover the assay plates with lids and incubate in a 37°C, 5% CO₂ incubator for 16–24 hours.

Note: The incubation period was optimized using blinatumomab. We recommend optimizing assay time (3–24 hours) with your biologic drug of interest.

8. Near the end of the assay incubation, equilibrate the Bio-Glo-NB™ TCK Luciferase Assay Buffer to ambient temperature.

5.H. Day 2: Preparing and Adding Bio-Glo-NB™ TCK Reagent

We recommend preparing the Bio-Glo-NB™ TCK Reagent immediately before use. Ensure that the Bio-Glo-NB™ TCK Luciferase Assay Buffer is equilibrated to room temperature (do not exceed 25°C) before reconstituting the reagent. Do not store the reconstituted reagent. Once reconstituted, the reagent will lose approximately 15% activity over 8 hours and approximately 60% activity over 24 hours at room temperature.

 **Note:** HiBiT TCK Bioassays are compatible only with Bio-Glo-NB™ TCK Reagent. **Do not** use Bio-Glo™ Reagent or Bio-Glo-NL™ Reagent.

1. Remove the Bio-Glo-NB™ TCK Luciferase Assay Substrate from –30°C to –10°C storage and mix by pipetting. Briefly centrifuge the tubes if the substrate has collected in the cap or on the sides of the tubes.
2. Remove the Bio-Glo-NB™ TCK LgBiT Protein from –30°C to –10°C storage and mix by pipetting. Briefly centrifuge the tubes if the solution has collected in the cap or on the sides of the tubes.
3. Transfer the desired amount of ambient temperature Bio-Glo-NB™ TCK Luciferase Assay Buffer to a 15ml or 50ml centrifuge tube.
4. Add Bio-Glo-NB™ TCK LgBiT Protein (1:100 dilution) and Bio-Glo-NB™ TCK Luciferase Assay Substrate (1:50 dilution) to the Bio-Glo-NB™ TCK Luciferase Assay Buffer. For example, if the experiment requires 10ml of reagent, add 100µl of Bio-Glo-NB™ TCK LgBiT Protein and 200µl of Bio-Glo-NB™ TCK Luciferase Assay Substrate to 10ml of Bio-Glo-NB™ TCK Luciferase Assay Buffer. Ten milliliters (10ml) of Bio-Glo-NB™ TCK Reagent is sufficient for 120 wells (two assay plates, using the inner 60 wells of each plate).
5. Prepare the Bio-Glo-NB™ TCK Reagent with 200µg/ml digitonin for detection of MR controls as follows:
 - a. Calculate the volume of Bio-Glo-NB™ TCK Reagent required for MR control wells.
 - b. Transfer the calculated volume of fully reconstituted Bio-Glo-NB™ TCK Reagent to a separate tube.
 - c. Add digitonin to a final concentration of 200µg/ml.
6. Remove assay plates from the incubator after the incubation period and equilibrate to room temperature for 10–15 minutes.
7. Using a manual multichannel pipette, add 75µl the Bio-Glo-NB™ TCK Reagent to wells treated with test or reference biologic drugs and SR controls, taking care to not create bubbles.

8. Add 75µl of the Bio-Glo-NB™ TCK Reagent with 200µg/ml digitonin (from Step 5) to MR control wells, taking care to not create bubbles.
9. Wait 10 minutes, then measure the luminescence in a GloMax® Discover System or a plate reader with glow-type luminescence reading capabilities. The luminescence intensity will decay gradually, with a signal half-life of approximately 120 minutes at room temperature.

Note: Varying the Bio-Glo-NB™ TCK Reagent incubation time will affect the raw RLU values but should not significantly change the EC₅₀ value and maximum fold induction.

5.I. Data Analysis

The HiBiT TCK Bioassay is a nonlytic bioassay that measures HiBiT released from dead or dying cells using Bio-Glo-NB™ TCK Reagent. The luminescent signal is expressed as RLU.

1. Calculate fold induction:

$$\text{Fold Induction} = \frac{\text{RLU}_{\text{Induced}}}{\text{RLU}_{\text{SR control}}}$$

2. Calculate percent specific lysis:

$$\text{Percent Specific Lysis} = \frac{\text{RLU}_{\text{Induced}} - \text{RLU}_{\text{SR control}}}{\text{RLU}_{\text{MR control}} - \text{RLU}_{\text{SR control}}} \times 100$$

Maximum release (MR) control wells contain Raji Cells (LDH-HiBiT) plus CD8+ T cells (no test or reference biologic drug added). These wells are treated with Bio-Glo-NB™ TCK Reagent supplemented with 200µg/ml digitonin to measure the total amount of HiBiT protein in each well.

Spontaneous release (SR) control wells contain Raji Cells (LDH-HiBiT) plus CD8+ T cells (no test or reference biologic drug added). These wells measure background killing in the absence of a test or reference biologic drug.

3. Graph data as RLU versus Log₁₀ [biologic drug], fold induction versus Log₁₀ [biologic drug] and percent specific lysis versus Log₁₀ [biologic drug]. Fit curves and determine the EC₅₀ value of the biologic drug response using appropriate curve fitting software (such as GraphPad Prism® software).

5.J. Using Knockout Controls

Knockout cells should be plated at the same density as the parental cell line and handled identically during bioassay setup. Users should continue to include test and reference antibodies on the same plate, using separate plates for parental and knockout cell lines. A biologic drug that targets an antigen remaining on the surface of knockout cells (e.g., BCMA) can be included as a positive control for both parental and knockout cell lines (Figure 6). Always calculate percent specific lysis when comparing between parental and knockout controls, as clones may vary in raw RLU.

6. CAR-T Assay Protocol

This procedure illustrates the use of the Raji (LDH-HiBiT) Bioassay to test CAR-T effector cells against a reference sample in a single cytotoxicity assay (Figure 13). Each test and reference CAR-T is assayed in triplicate, in an eight-point effector-to-target (E:T) ratio dilution series, in a single 96-well assay plate using the inner 60 wells. Other experimental and plate layouts are possible but may require further optimization. See Section 6.I for recommendations on including knockout cell line controls.

Note: When preparing test and reference CAR-T dilutions, choose an appropriate starting concentration and dilution scheme to achieve a complete dose-response curve with proper upper and lower asymptotes and sufficient points on the slope. For reference, we use a 100:1 E:T ratio as a starting concentration (1X) and twofold serial dilutions when testing CAR-T cells. Appropriate dilution schemes may vary significantly depending on the potency of your CAR-T cells.

Maximum release (MR) and spontaneous release (SR) controls should be included in your experiment. The MR control utilizes digitonin to permeabilize target cells, providing the maximal luminescent signal. The SR control is a measure of background cell death in the absence of CAR-T cells.

Materials to Be Supplied By the User

(Composition of Buffers and Solutions is provided in Section 9.A.)

Reagents

- user-provided test and reference CAR-T effector cells recognizing antigen expressed on Raji Cells (LDH-HiBiT) (CD19, CD20, CD22, etc.)
- digitonin (e.g., Cat.# G9441)

Supplies and Equipment

- solid-white, U-bottom 96-well assay plates (e.g., Corning® Cat.# 3355)
- sterile clear V-bottom 96-well plate with lid (e.g., Costar® Cat.# 3896) for preparing CAR-T cell dilutions
- pipettes (single-channel and 12-channel; for best results use both manual and electronic pipettes as needed)
- sterile 15ml and 50ml conical tubes
- sterile reagent reservoirs (e.g., Corning®/Costar® Cat.# 4870)
- 37°C, 5% CO₂ incubator
- 37°C water bath
- T25 cell culture flask (e.g., Corning Cat #430639)
- sensitive plate reader with luminescence measuring capability or luminometer (e.g., GloMax® Discover System or equivalent system)

Notes:

- a. Solid-white flat-bottom assay plates can also be used for CAR-T assays. However, with flat-bottom plates we generally observe a right-shifted dose-response curve compared to U-bottom assay plates. Optimal plate type will need to be evaluated for each CAR-T product.
- b. Corning® Cat.# 3355 plates are provided without lids and are marked as nonsterile. Use a lid from another Corning® plate type (e.g., Cat.# 3917). Although marked as nonsterile, we had no issues with sterility when using these plates while developing this product.




Figure 13. Schematic protocol for the Raji (LDH-HiBiT) TCK Bioassay in a CAR-T assay application.

6.A. Preparing Assay Reagents

Heat-Inactivated (HI) FBS: Thaw the provided FBS overnight at +2°C to +10°C or in a 37°C water bath on the day of use. Heat inactivate the FBS for 30 minutes at 56°C. Mix well and adjust to 37°C before use.

Raji Cells (LDH-HiBiT) Recovery Medium: On the day before the assay, prepare 15ml of Raji Cells (LDH-HiBiT) recovery medium (see Section 9.A). Thaw both vials of FBS overnight at +2°C to +10°C or in a 37°C water bath and heat inactivate as noted above. To prepare Raji Cells (LDH-HiBiT) recovery medium, add 1.5ml of FBS to 13.5ml of 37°C RPMI 1640 Medium. For reference, 15ml of Raji Cells (LDH-HiBiT) recovery medium is sufficient to thaw and plate one vial of Raji cells (LDH-HiBiT). If multiple vials will be thawed, then scale the amount of Raji Cells (LDH-HiBiT) recovery medium appropriately. Store the remaining HI-FBS at +2°C to +10°C for use in preparing the assay buffer on the day of the assay.

CAR-T Assay Buffer: On the day of the assay, prepare an appropriate amount of CAR-T assay buffer (90% RPMI 1640 with 10% HI-FBS). Mix well and warm to 37°C before use.

Note: The recommended assay buffer contains 10% HI-FBS. This concentration of FBS works well for the CAR-T that we have tested. If you experience assay performance issues when using this assay buffer, we recommend testing different serum concentrations in the range of 0.5–10%.

Bio-Glo-NB™ TCK Reagent: For reference, 10ml of Bio-Glo-NB™ TCK Reagent is sufficient to assay 120 wells in a 96-well assay format. Store the Bio-Glo-NB™ TCK Luciferase Assay Substrate and Bio-Glo-NB™ TCK LgBiT Protein at –30°C to –10°C. Thaw the Bio-Glo-NB™ TCK Luciferase Assay Buffer at room temperature (do not exceed 25°C) near the end of the assay induction period. We recommend preparing the reconstituted Bio-Glo-NB™ TCK Reagent immediately before use.



Note: HiBiT TCK Bioassays are compatible only with Bio-Glo-NB™ TCK Reagent. **Do not** use Bio-Glo™ Reagent or Bio-Glo-NL™ Reagent.

Test and Reference Samples: Using assay buffer as the diluent, prepare stock starting dilutions (dilu1, 2X final concentration of CAR-T cells for the highest E:T ratio) of each CAR-T before making serial dilutions.

6.B. Plate Layout Design

For the protocol described here, use the plate layout illustrated in Figure 14 as a guide. The protocol describes serial replicate dilutions ($n = 3$) of test and reference CAR-T cell preparations to generate two eight-point dose-response curves for each plate.

Recommended Plate Layout Design													
	1	2	3	4	5	6	7	8	9	10	11	12	
A	B	B	B	B	B	B	B	B	B	B	B	B	Assay Buffer (B)
B	B	SR	dilu8	dilu7	dilu6	dilu5	dilu4	dilu3	dilu2	dilu1	MR	B	Reference CAR-T
C	B	SR	dilu8	dilu7	dilu6	dilu5	dilu4	dilu3	dilu2	dilu1	MR	B	Test CAR-T
D	B	SR	dilu8	dilu7	dilu6	dilu5	dilu4	dilu3	dilu2	dilu1	MR	B	Reference CAR-T
E	B	SR	dilu8	dilu7	dilu6	dilu5	dilu4	dilu3	dilu2	dilu1	MR	B	Test CAR-T
F	B	SR	dilu8	dilu7	dilu6	dilu5	dilu4	dilu3	dilu2	dilu1	MR	B	Reference CAR-T
G	B	SR	dilu8	dilu7	dilu6	dilu5	dilu4	dilu3	dilu2	dilu1	MR	B	Test CAR-T
H	B	B	B	B	B	B	B	B	B	B	B	B	Assay Buffer (B)

Figure 14. Example plate layout showing nonclustered sample locations of test and reference CAR-T dilution series.

Maximum release (MR) and spontaneous release (SR) control wells contain Raji Cells (LDH-HiBiT) alone (no effector cells). MR control wells receive Bio-Glo-NB™ TCK Reagent containing 200µg/ml digitonin. Wells containing assay buffer alone are denoted by “B”.

6.C. Thawing Raji Cells (LDH-HiBiT)

Thaw-and-use Raji Cells (LDH-HiBiT) are in a sensitive state immediately post-thaw and care should be taken to follow the thawing and plating procedures exactly as described. Handle the cells gently and do not overmix the cells or overwarm the reagents. Do not vortex the cells or centrifuge at speeds $>150 \times g$. It is important to keep pipette tips away from tube surfaces during delivery or mixing steps to minimize mechanical stress on the cells. We recommend that you thaw and dilute a maximum of two vials of thaw-and-use cells at any one time.

6.C. Thawing Raji Cells (LDH-HiBiT) (continued)

! **Note:** Perform the following steps in a sterile cell culture hood.

1. Add 9.5ml of Raji Cells (LDH-HiBiT) recovery medium to a 15ml conical tube and warm to 37°C.
2. Remove one cryovial of Raji Cells (LDH-HiBiT) from storage at -140°C and thaw in a 37°C water bath with gentle agitation (do not invert) until just thawed (typically 2–3 minutes).
3. Using a 2ml serological pipette or a 1ml wide-bore pipette tip, transfer 1ml of prewarmed Raji Cells (LDH-HiBiT) recovery medium, one drop at a time, to the cryovial. Using the same pipette, gently mix the contents of the cryovial by slowly pipetting approximately 1ml, once. Gently transfer the contents of the cryovial to the 15ml conical tube containing prewarmed Raji Cells (LDH-HiBiT) recovery medium.
4. Pellet the cells at 150 × g for 5 minutes.
5. Aspirate the supernatant, leaving 100–300µl of Raji Cells (LDH-HiBiT) recovery medium above the cell pellet.
6. Resuspend the cell pellet by gently flicking the tube bottom 2–3 times. Add 5ml of prewarmed Raji Cells (LDH-HiBiT) recovery medium to the side of the tube. Do not pipet Raji Cells (LDH-HiBiT) recovery medium directly onto the cell pellet. Pipet once to mix.
7. Transfer the cell suspension to a T25 cell culture flask and place the flask horizontally in a humidified 37°C, 5% CO₂ incubator.
8. Incubate for approximately 16–24 hours before assay setup.

6.D. Day 2: Preparing Raji Cells (LDH-HiBiT)

! **Note:** Perform the following steps in a sterile cell culture hood.

1. Gently mix and harvest the Raji Cells (LDH-HiBiT) that were cultured overnight.
2. Transfer into a 50ml centrifuge tube.
3. Centrifuge cells at 150 × g for 5 minutes.
4. Wash cells by removing supernatant, and gently resuspending in 10ml of prewarmed assay buffer.
5. Centrifuge cells at 150 × g for 5 minutes, then remove and discard supernatant.
6. Gently resuspend the cell pellet in 16ml of prewarmed assay buffer.
7. Incubate the Raji Cells (LDH-HiBiT) at ambient temperature until Section 6.F.

6.E. Day 2: Preparing 2X Stock Solutions of CAR-T Cells

1. Thaw or harvest CAR-T cells according to the protocol established in your lab.
2. Pellet cells according to your protocol and resuspend in CAR-T assay buffer at a concentration of 6.25 × 10⁶ viable cells/ml, a 2X stock solution for a maximum E:T ratio of 50:1.

3. In a clear V-bottom 96 well plate or 12-well dilution reservoir, perform seven twofold serial dilutions of CAR-T cells using CAR-T assay buffer as diluent. Prepare enough of each cell suspension for 40µl per well of each cell density tested in triplicate.
4. Proceed immediately to Section 6.F.

6.F. Day 2: Adding Target Cells and CAR-T Cells to Assay Plates


1. Gently resuspend and transfer the Raji Cells (LDH-HiBiT) suspension to a sterile reagent reservoir.
2. Using a multichannel pipette, immediately dispense 40µl of the cell suspension to each of wells B2–G11 of a 96-well white U-bottom assay plate. The final cell concentration should be 5×10^3 cells/well.
3. Using a multichannel pipette, add 40µl of the appropriate CAR-T dilution to wells B3–G10 of the assay plates according to the plate layout in Figure 14. Add 40µl of CAR-T assay buffer to MR and SR control wells (B2–G2 and B11–G11, respectively).

Note: Digitonin will be added at the time of luminescence measurement (Section 6.G).

4. Add 80µl of CAR-T assay buffer to each of the outside wells of the assay plates.
5. Cover the assay plates with lids and incubate in a 37°C, 5% CO₂ incubator for 4–72 hours. We recommend testing several assay times as part of assay optimization for the CAR-T of interest.
6. Near the end of the assay incubation, warm Bio-Glo-NB™ TCK Luciferase Assay Buffer to ambient temperature.

6.G. Day 2: Preparing and Adding Bio-Glo-NB™ TCK Reagent

We recommend preparing the Bio-Glo-NB™ TCK Luciferase Reagent immediately before use. Ensure that the Bio-Glo-NB™ TCK Luciferase Assay Buffer is equilibrated to room temperature (do not exceed 25°C) before reconstituting the reagent. Do not store the reconstituted reagent. Once reconstituted, the reagent will lose approximately 15% activity over 8 hours and approximately 60% activity over 24 hours at room temperature.

 **Note:** HiBiT TCK Bioassays are compatible only with Bio-Glo-NB™ TCK Reagent. **Do not** use Bio-Glo™ Reagent or Bio-Glo-NL™ Reagent.

1. Remove the Bio-Glo-NB™ TCK Luciferase Assay Substrate from –30°C to –10°C storage and mix by pipetting. Briefly centrifuge the tubes if substrate has collected in the cap or on the sides of the tubes.
2. Remove the Bio-Glo-NB™ TCK LgBiT Protein from –30°C to –10°C storage and mix by pipetting. Briefly centrifuge the tubes if the solution has collected in the cap or on the sides of the tubes.
3. Transfer the desired amount of ambient temperature Bio-Glo-NB™ TCK Luciferase Assay Buffer to a 15ml or 50ml centrifuge tube.
4. Add Bio-Glo-NB™ TCK LgBiT Protein (1:100 dilution) and Bio-Glo-NB™ TCK Luciferase Assay Substrate (1:50 dilution) to the Bio-Glo-NB™ TCK Luciferase Assay Buffer. For example, if the experiment requires 10ml of reagent, add 100µl of Bio-Glo-NB™ TCK LgBiT Protein and 200µl of Bio-Glo-NB™ TCK Luciferase Assay Substrate to 10ml of Bio-Glo-NB™ TCK Luciferase Assay Buffer. Ten milliliters (10ml) of the Bio-Glo-NB™ TCK Reagent is sufficient for 120 wells (two assay plates, using the inner 60 wells of each plate).

6.G. Day 2: Preparing and Adding Bio-Glo-NB™ TCK Reagent (continued)

5. Prepare Bio-Glo-NB™ TCK Reagent with 200µg/ml digitonin as a detection reagent for MR controls:
 - a. Calculate the volume of Bio-Glo-NB™ TCK Reagent required for the MR control wells.
 - b. Transfer the calculated volume of Bio-Glo-NB™ TCK Reagent to a separate tube.
 - c. Add digitonin to a final concentration of 200µg/ml.
 6. Remove assay plates from the incubator after the incubation period and equilibrate to room temperature for 10–15 minutes.
 7. Using a manual multichannel pipette, add 80µl of Bio-Glo-NB™ TCK Reagent to the wells treated with test or reference CAR-T cells and SR controls, taking care not to create bubbles.
 8. Add 80µl of Bio-Glo-NB™ TCK Reagent plus 200µg/ml of digitonin (from Step 5) to MR control wells.
 9. Wait 10 minutes, then measure the luminescence in a GloMax® Discover System or a plate reader with glow-type luminescence reading capabilities. The luminescence intensity will decay gradually, with a signal half-life of approximately 1–2 hours at room temperature.
- Note:** Varying the Bio-Glo-NB™ TCK Reagent™ incubation time will affect the raw relative light unit (RLU) values but should not significantly change the EC₅₀ value and maximum fold induction.

6.H. Data Analysis

The HiBiT TCK Bioassay is a nonlytic bioassay that measures HiBiT released from dead or dying cells using Bio-Glo-NB™ TCK Reagent. Luminescence signal is expressed as RLU.

1. Calculate fold induction:

$$\text{Fold Induction} = \frac{\text{RLU}_{\text{Induced}}}{\text{RLU}_{\text{SR control}}}$$

2. Calculate percent specific lysis:

$$\text{Percent Specific Lysis} = \frac{\text{RLU}_{\text{Induced}} - \text{RLU}_{\text{SR control}}}{\text{RLU}_{\text{MR control}} - \text{RLU}_{\text{SR control}}} \times 100$$

Maximum release (MR) control wells contain Raji Cells (LDH-HiBiT) alone (no CAR-T cells added). These wells are treated with Bio-Glo-NB™ TCK Reagent supplemented with 200µg/ml digitonin to measure the total amount of HiBiT protein in each well.

Spontaneous release (SR) control wells contain Raji Cells (LDH-HiBiT) alone (no CAR-T cells added). These wells measure background target cell death in the absence of CAR-T cells.

3. Graph data as RLU versus Log₁₀(E:T ratio), fold induction versus Log₁₀(E:T ratio) and percent specific lysis versus Log₁₀(E:T ratio). Fit curves and determine the EC₅₀ value of the CAR-T response using appropriate curve fitting software (such as GraphPad Prism® software).

6.1. Using Knockout Controls

Knockout cells should be plated at the same density as the parental cell line and handled identically during bioassay setup. Users should continue to include test and reference preparations of CAR-T cells on the same plate, using separate plates for parental and knockout cell lines. A CAR-T preparation that targets an antigen remaining on the surface of knockout cells (e.g., CD22 or BCMA) can be included as a positive control for both parental and knockout cell lines. Always calculate percent specific lysis when comparing between parental and knockout controls, as clones may vary in raw RLU values.

7. Troubleshooting

For questions not addressed here, please contact your local Promega Branch Office or Distributor. Contact information available at: www.promega.com. Email: techserv@promega.com

Symptoms	Causes and Comments
Low luminescence measurements (RLU readout)	<p>Ensure that you are using Bio-Glo-NB™ TCK Reagent. Raji (LDH-HiBiT) TCK Bioassays are not compatible with Bio-Glo™ or Bio-Glo-NL™ Reagents.</p> <hr/> <p>Choose an instrument designed for plate-reading luminescence detection. Instruments designed primarily for fluorescence detection are not recommended. Luminometers measure and report luminescence as relative values, and actual RLU numbers will vary between instruments.</p> <hr/> <p>Some models of luminometers with low sensitivity should be avoided. If using a reader with an adjustable gain, we recommend a high-gain setting.</p> <hr/> <p>Insufficient cells per well can lead to low RLU. Handle and plate cells according to the instructions to ensure a sufficient number of viable cells per well.</p> <hr/> <p>Low activity of Bio-Glo-NB™ TCK Reagent leads to low RLU. Store and handle Bio-Glo-NB™ Reagent according to the instructions. For best results, prepare immediately before use.</p>
Weak assay response (low fold induction)	<p>Ensure that you are using Bio-Glo-NB™ TCK Reagent. Raji (LDH-HiBiT) TCK Bioassays are not compatible with Bio-Glo™ or Bio-Glo-NL™ Reagents.</p> <hr/> <p>Optimize the concentration range of your test sample(s) to achieve a full dose response with complete upper and lower asymptotes. The EC₅₀ value obtained with the HiBiT TCK Bioassay may vary from the EC₅₀ value obtained using other methods.</p>

7. Troubleshooting (continued)

Symptoms	Causes and Comments
Weak assay response (low fold induction, continued)	<p>The assay is sensitive to the concentration of FBS in assay buffer. Optimize the FBS concentration from 0.5–10% in assay buffer if assay performance is not ideal. Always use heat-inactivated FBS.</p> <p>Optimize the assay incubation time within a range of 3–24 hours (or longer for the CAR-T assay).</p> <p>Performance of PBMC ADCC assays is highly donor-dependent. Target cells associated with Raji (LDH-HiBiT) TCK Bioassays have been optimized for use with Promega Human PBMC ADCC-Qualified.</p> <p>Always use white, round-bottom plates for PBMC ADCC assays. In most cases, TDCC and CAR-T assays can be performed in round- or flat-bottom plates, although the performance characteristics of your biologic drug or engineered effector cells may vary with plate type.</p>
High spontaneous release	<p>If spontaneous release RLU is less than 100-fold above plate reader background RLU, subtract plate background RLU from all samples before calculating fold induction.</p> <p>Insufficient washing of target cells prior to assay setup leads to high spontaneous release values due to carryover of free HiBiT from culture. Always wash target cells thoroughly, according to the protocols listed in this technical manual.</p> <p>Low cell viability causes high spontaneous release values because dead and dying cells release HiBiT into the medium. Target cells should be >95% viable at the time of assay setup. Handle cells carefully and according to the instructions in this protocol. Centrifuge at low speeds and resuspend cell pellets gently.</p>
Specific lysis greater than 100%	<p>In assays longer than 6 hours, Raji Cells (LDH-HiBiT) may divide before being lysed by effector cells, leading to an increase in the total HiBiT present in the well. In short (same-day) assays, add digitonin to maximum release wells at the time of assay setup. For overnight assays (or longer) add digitonin when adding detection reagent.</p> <p>Digitonin at concentrations greater than 200µg/ml and alternative detergents can interfere with NanoBiT® complementation and reduce luminescence in maximum release wells. Always use digitonin at the recommended final concentration of 100µg/ml.</p>

Symptoms

Variability in assay performance

Causes and Comments

Inappropriate handling of target cells, including excessive centrifugation times and speeds, can cause low assay performance. Centrifuge exactly according to instructions.

Inappropriate cell counting methods can lead to variation in cell numbers in culture and assays and cause high assay variation. Ensure consistent and accurate cell counting methods.

RLU differ slightly between parental and knockout cell lines

Use percent specific lysis to compare the functional response between the parental clone and one or more knockout clones if the small difference in light output the different cell lines is a concern.

8. References

1. Barnhart, B.C. and Quigley, M. (2017) Role of Fc-FcγR interactions in the antitumor activity of therapeutic antibodies. *Immunol. Cell Biol.* **95**, 340–6.
2. van de Donk, N.W.C.J. and Zweegman, S. (2023) T-cell-engaging bispecific antibodies in cancer. *Lancet* **402**, 142–58.
3. Haslauer, T. *et. al.* (2021) CAR T-cell therapy in hematological malignancies. *Int. J. Mol. Sci.* **22**, 8996.
4. Schwinn, M.K. *et. al.* (2018) CRISPR-mediated tagging of endogenous proteins with a luminescent peptide. *ACS Chem. Biol.* **13**, 467–74.

9. Appendix

9.A. Composition of Buffers and Solutions

recovery medium

90% RPMI 1640 with L-glutamine and HEPES
10% HI-FBS

CAR-T assay buffer

90% RPMI 1640 with L-glutamine and HEPES
10% HI-FBS

PBMC ADCC assay buffer

95% RPMI 1640 with L-glutamine and HEPES
5% HI-FBS
55µM β-mercaptoethanol
5ng/ml IL-2

PBMC culture medium

90% RPMI 1640 with L-glutamine and HEPES
10% HI-FBS
1X sodium pyruvate
55µM β-mercaptoethanol
5ng/ml IL-2

TDCC assay buffer

90% RPMI 1640 with L-glutamine and HEPES
10% HI-FBS

9.B. Related Products

Fc Effector Bioassays

Product	Size	Cat.#
ADCC Reporter Bioassay, Complete Kit (Raji)*	1 each	G7015
ADCC Reporter Bioassay, Core Kit*	1 each	G7010
ADCC Reporter Bioassay, F Variant, Core Kit**	1 each	G9790
ADCC Reporter Bioassay, Target Kit (Raji)*	1 each	G7016
FcγRIIIa-H ADCP Reporter Bioassay, Complete Kit**	1 each	G9901
FcγRIIIa-H ADCP Reporter Bioassay, Core Kit**	1 each	G9991
Mouse FcγRIV ADCC Bioassay, Complete Kit**	1 each	M1201
Mouse FcγRIV ADCC Bioassay, Core Kit**	1 each	M1211
Membrane TNFα Target Cells**	1 each	J3331
Membrane RANKL Target Cells**	1 each	J3381

*For Research Use Only. Not for use in diagnostic procedures.

**Not for Medical Diagnostic Use.

Additional kit formats are available.

9.B. Related Products (continued)

Fc Effector Immunoassay

Product	Size	Cat.#
Lumit™ FcRn Binding Immunoassay	100 assays	W1151

Not for Medical Diagnostic Use. Additional kit formats and sizes are available.

Immune Checkpoint Bioassays

Product	Size	Cat.#
4-1BB Bioassay	1 each	JA2351
CD28 Bioassay	1 each	JA6701
CD28 Blockade Bioassay	1 each	JA6101
CD40 Bioassay	1 each	JA2151
CTLA-4 Blockade Bioassay	1 each	JA3001
GITR Bioassay	1 each	JA2291
ICOS Bioassay	1 each	JA6801
ICOS Blockade Bioassay	1 each	JA6001
LAG-3/MHCII Blockade Bioassay	1 each	JA1111
OX40 Bioassay	1 each	JA2191
PD-1/PD-L1 Blockade Bioassay	1 each	J1250
PD-1+TIGIT Combination Bioassay	1 each	J2211
PD-L1 Negative Cells	1 each	J1191
TIGIT/CD155 Blockade Bioassay	1 each	J2201

Not for Medical Diagnostic Use. Additional kit formats and sizes are available.

T Cell Activation Bioassays

Product	Size	Cat.#
T Cell Activation Bioassay (IL-2)	1 each	J1651
T Cell Activation Bioassay (NFAT)	1 each	J1621
T Cell Activation Bioassay (TCRαβ-KO, CD4+)	1 each	GA1172
T Cell Activation Bioassay (TCRαβ-KO, CD8+)	1 each	GA1162
T Cell Activation Bioassay (TCRαβ-KO, CD4+, CD8+)	1 each	GA1182

Not for Medical Diagnostic Use. Additional kit formats and sizes are available.

Cytokine and Growth Factor Bioassays

Product	Size	Cat.#
IL-2 Bioassay	1 each	JA2201
IL-6 Bioassay	1 each	JA2501
IL-12 Bioassay	1 each	JA2601
IL-15 Bioassay	1 each	JA2011
IL-23 Bioassay	1 each	JA2511
RANKL Bioassay	1 each	JA2701
VEGF Bioassay	1 each	GA2001

Not for Medical Diagnostic Use. Additional kit formats and sizes are available.

Macrophage-Directed Bioassays

Product	Size	Cat.#
SIRP α /CD47 Blockade Bioassay	1 each	JA6011
SIRP α /CD47 Blockade Bioassay, Fc-Dependent	1 each	JA4801
TLR Bioassay	1 each	JA9011
ADCP Reporter Bioassay (THP-1)	1 each	JA9411

Not for Medical Diagnostic Use. Additional kit formats are available.

HiBiT Target Cell Killing Bioassays

Product	Size	Cat.#
Raji (HT-HiBiT) TCK Bioassay	1 each	JA1211
Raji (LDH-HiBiT) TCK Bioassay	1 each	JA1311
Ramos (HiBiT) TCK Bioassay	1 each	JA1411
H929 (HiBiT) TCK Bioassay	1 each	JA1511

Not for Medical Diagnostic Use. Additional kit formats are available.

Control Antibodies and Proteins

Product	Size	Cat. #
Control Ab, Anti-4-1BB	50µg	K1161
Control Ab, Anti-CD20	5µg	GA1130
Control Ab, Anti-CD40	50µg	K1181
Control Ab, Anti-CTLA-4	100µg	JA1020
Control Ab, Anti-LAG-3	100µg	K1150
Control Ab, Anti-OX40	50µg	K1191
Control Ab, Anti-PD-1	100µg	J1201
Control Ab, Anti-SIRPα	50µg	K1251
Control Ab, Anti-TIGIT	100µg	J2051
Control Ab, Anti-TIM-3	100µg	K1210
Recombinant VEGF ligand	10µg	J2371

Detection Reagents

Product	Size	Cat. #
Bio-Glo™ Luciferase Assay System	10ml	G7941
Bio-Glo-NL™ Luciferase Assay System	10ml	J3081
Bio-Glo-NB™ TCK Luciferase Assay System	10ml	JB1001

Not for Medical Diagnostic Use. Additional sizes are available.

Detection Instruments

Product	Size	Cat. #
GloMax® Navigator System	1 each	GM2000
GloMax® Discover System	1 each	GM3000
GloMax® Explorer System	1 each	GM3500

For Research Use Only. Not for use in diagnostic procedures.

Note: Additional Fc Effector, Immune Checkpoint, T Cell Activation, Cytokine, Macrophage, Primary Cell and Target Cell Killing Bioassays are available. To view and order Promega Bioassay products visit:

www.promega.com/products/reporter-bioassays/ or email: EarlyAccess@promega.com. For information on custom bioassay development and services visit the Promega Tailored R&D Solutions website:

www.promega.com/custom-solutions/tailored-solutions/



©NOT FOR MEDICAL DIAGNOSTIC USE. FOR IN VITRO USE ONLY. BY USE OF THIS PRODUCT, RECIPIENT AGREES TO BE BOUND BY THE TERMS OF THIS LIMITED USE LABEL LICENSE. If recipient is not willing to accept the terms of this label license, and the product is unused, Promega will accept return of the unused product and provide recipient with a full refund.

This product may not be further sold or transferred by the recipient and may be used only by the recipient, and then only for research use, which includes but is not limited to drug discovery and development, quality assurance testing, and product release assays. No other commercial use is allowed. "Commercial use" means any and all uses of this product by recipient for monetary or other consideration, including providing a service, information or data to unaffiliated third parties, and resale of this product for any use. Recipient has no right to modify, derivatize, genetically engineer or otherwise create variations of the product, including but not limited to the nucleotide sequence encoding the HiBiT peptide, for which the sequence has been stably engineered within the cells. Recipient may propagate and store the cells for long-term use. In addition, recipient must use Bio-Glo-NB™ TCK or Bio-Glo-NB™ Lytic Luciferase Assay Systems purchased from Promega for all luminescence assays using this product or contact Promega to obtain a license for use of this product with reagents other than Promega's.

PROMEGA MAKES NO REPRESENTATIONS OR WARRANTIES OF ANY KIND, EITHER EXPRESSED OR IMPLIED, INCLUDING FOR MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE WITH REGARDS TO THIS PRODUCT. The terms of this label license shall be governed under the laws of the State of Wisconsin, USA.

©By purchase of this material, Promega Corporation ("Promega") conveys to purchaser (the "Limited Licensee") the non-transferable right to use the material purchased from Promega and derivatives made by Limited Licensee solely for research conducted by Limited Licensee in accordance with the following requirements: (i) Limited Licensee shall not sell or otherwise transfer material or derivatives to any other person or entity, or use material or derivatives to perform services for the benefit of any other person or entity, except to transfer the material or derivatives to a Limited Licensee affiliate provided such affiliate agrees to be bound by the terms of this Limited License; (ii) Limited Licensee shall use material and derivatives only for its internal research within the Customer Field, which may include internal research within the Customer Field in connection with product research, and not for any Commercial Purposes; (iii) Limited Licensee shall use material and derivatives in compliance with all applicable laws and regulations, including without limitation applicable human health and animal welfare laws and regulations; (iv) Institutions shall provide no warranties of any kind to Limited Licensee (statutory or implied) concerning the Intellectual Property or material or derivatives, including without limitation, as to product quality, condition, description, merchantability, fitness for a particular purpose, noninfringement of intellectual property rights or the absence of latent or other defects, and all such warranties are hereby expressly disclaimed; (v) Institutions shall expressly disclaim any warranty regarding results obtained through the use of material or derivatives, including without limitation any claim of inaccurate, invalid or incomplete results; (vi) Institutions and their directors, trustees, officers, employees, agents, faculty, affiliated investigators and students shall have no liability to Limited Licensee, including without limitation, for any loss of use or profits, business interruption or any consequential, incidental, special or other indirect damages of any kind, regardless of how caused and regardless of whether an action in contract, tort, strict product liability or otherwise; (vii) Limited Licensee shall indemnify, defend and hold harmless the Indemnitees and HHMI Indemnitees against any liability, damage, loss or expense (including without limitation reasonable attorneys' fees and expenses) incurred by or imposed upon any of the Indemnitees or HHMI Indemnitees, as applicable, in connection with any claims, suits, investigations, actions, demands or judgments arising out of or related to the exercise of any rights granted to Limited Licensee under this limited license or any breach of limited license by such Limited Licensee, provided that, to the extent the foregoing is not permitted by law, Limited Licensee agrees, to the extent permitted by law, that it, and not the Indemnitees or HHMI Indemnitees, as applicable, shall be responsible for any liability, damage, loss or expense arising out of or related to the exercise of any rights granted to Limited Licensee under this limited license or any breach of this limited license by Limited Licensee; and (viii) material and derivatives and their use may be the subject of one or more issued patents or pending patent applications owned by one or more Institutions, and the purchase of material does not convey a license under any claims in the foregoing patents or patent applications directed to the material or derivatives or use, production or commercialization thereof, except as expressly set forth in this limited license. Nothing in this limited license shall be construed to confer any rights upon Limited Licensee by implication, estoppel or otherwise as to any technology or patent rights of Broad or any other entity other than the Intellectual Property. In addition, nothing in this limited license shall be construed to confer upon Limited Licensee or any third party any rights under or to the Intellectual Property outside of the Customer Field.

Commercial Purposes means (a) the practice, performance or provision of any methods, process or service; or (b) the manufacture, production, sale, use, distribution, disposition or importing of any product, in each case (a) or (b) for consideration of any kind, for the purpose of sale or commercial exploitation, or on any other commercial basis.

Customer Field means use as a research tool for research purposes and use for quality control/quality assurance testing and product lot release assays, provided, however, that notwithstanding the foregoing, the Customer Field shall expressly exclude:

1. Any human or clinical use, including, without limitation, any administration into humans or any diagnostic or prognostic use;
2. Any human germline modification, including modifying the DNA of human embryos or human reproductive cells;
3. Any in vivo veterinary or livestock use;
4. The development, manufacture (including any bioproduction), distribution, importation, exportation, transportation, sale, offer for sale, marketing, promotion or other exploitation or use of the Intellectual Property or the material or derivatives as a therapeutic or diagnostic for humans or animals;
5. Products that provide nutritional benefits and are regulated by a regulatory authority as a drug or biologic pursuant to Section 505 of the Federal Food, Drug, and Cosmetic Act of 1938, as amended, Section 351 of the Public Health Service Act of 1944, as amended, or any successor laws, or equivalent laws or regulations in jurisdiction outside the United States;

6. Any agricultural use, including but not limited to, the use or application in the cultivation, growth, manufacture, exportation, or production of any tobacco product; and
7. Any use or application relating to gene drive.

HHMI Indemnitees means Howard Hughes Medical Institute and its trustees, officers, employees and agents.

Indemnitees means each Institution, their affiliates, and their current and former trustees, directors, officers, faculty, affiliated investigators, students, employees, medical and professional staff and agents and their respective successors, heirs and assigns.

Institutions means each of The Broad Institute, Inc., Massachusetts Institute of Technology, Harvard College, University of Iowa Research Foundation, University of Tokyo, The Rockefeller University, New York University, New York Genome Center, and Whitehead Institute of Biomedical Research collectively.

Intellectual Property means (a) any patents and patent applications listed on Appendix A; (b) any substitutions, divisions, continuations, reissues, renewals, re-examinations or extensions of any patents or patent applications of (a), any continuations-in-part of any patents or patent applications of (a) to the extent the claims are directed to subject matter specifically described in a patent or patent application set forth on Appendix A, and any application claiming priority to a provisional application set forth on Appendix A to the extent the claims are directed to subject matter specifically described in a provisional patent application set forth on Appendix A; (c) any letters patents and/or the equivalent issuing therefrom in any jurisdiction based on, and only to the extent of, any of the foregoing; and (d) any foreign or international equivalents of any of the foregoing, including without limitation any applications filed under the Patent Cooperation Treaty.

Appendix A, Broad IP List

BROAD TECHID	APPLICANTS	EXEMPLARY FAMILY SERIAL NO
BI-2011/008	Broad, MIT, Harvard and Rockefeller	PCT/US2013/074611; WO2014093595
BI-2011/008B	Broad, MIT and Harvard	14/259,420
BI-2011/008A	Broad and MIT	PCT/US2013/074743; WO2014093661
BI-2011/008C	Broad, MIT, and Harvard	PCT/US2013/074790; WO2014093694
BI-2011/020	Broad, MIT and Harvard	PCT/US2013/051418; WO2014018423
BI-2012/084A	Broad and MIT	PCT/US2013/074825; WO2014093718
BI-2012/084B	Broad, MIT and Harvard	PCT/US2013/074812; WO2014093709
BI-2013/003	Broad, MIT and Harvard	PCT/US2013/074667; WO2014093622
BI-2013/003D	Broad, MIT and Harvard	PCT/US2014/041803; WO2014201425
BI-2013/004E	Broad, MIT and Harvard	PCT/US2013/074691; WO2014093635
BI-2013/004F	Broad, MIT and Harvard	PCT/US2013/074736; WO2014093655
BI-2013/004G	Broad, MIT and Harvard	PCT/US2013/074819; WO2014093712
BI-2013/007	Broad, MIT and Harvard	14/855,046; US20160068822
BI-2013/066	Broad, MIT and Harvard	PCT/US2014/041800; WO2014204724
BI-2013/073	Broad and MIT	PCT/US2014/041806; WO2014204727
BI-2013/085	Broad, MIT and Whitehead	15/141,348; US20160251648
BI-2013/087J	Broad, Editas*, Iowa and MIT	PCT/US2014/064663; WO2015070083
BI-2013/087M	Broad, Iowa and MIT	PCT/US2014/069902; WO2015089354
BI-2013/087V	Broad, Iowa and MIT	PCT/US2014/069897; WO2015089351
BI-2013/093	Broad, MIT and Tokyo	15/171,141; US20160340660
BI-2013/094	Broad, MIT and Rockefeller	PCT/US2014/070135; WO2015089465
BI-2013/098	Broad and MIT	PCT/US2014/070068; WO2015089427

BI-2013/099	Broad, MIT and Harvard	PCT/US2014/041804; WO2014204726
BI-2013/101	Broad and MIT	PCT/US2014/070127; WO2015089462
BI-2013/103	Broad and MIT	PCT/US2014/041808; WO201404728
BI-2013/105	Broad, MIT and Harvard	PCT/US2014/041809; WO2014204729
BI-2013/107	Broad and MIT	PCT/US2014/070057; WO2015089419
BI-2013/112	Broad and MIT	PCT/US2013/074800; WO2014093701
BI-2013/113	Broad, MIT and Harvard	PCT/US2014/070152; WO2015089473
BI-2014/005	Broad, MIT, Harvard and Tokyo	PCT/US2014/070175; WO2015089486
BI-2014/061	Broad, MIT and Harvard	PCT/US2015/045504; WO2016028682
BI-2014/069	Broad and MIT	15/467,888; US20180010134
BI-2014/071	Broad and MIT	15/349,603; US20170107536
BI-2014/072	Broad and MIT	15/467,949; US20180044662
BI-2014/084	Broad, MIT and Harvard	15/469,081; US20180057810
BI-2014/097	Broad, MIT and Harvard	PCT/US2015/067177; WO2016106244
BI-2014/100	Broad, MIT and Harvard	PCT/US2015/065385; WO2016094867
BI-2014/101	Broad, MIT and Harvard	15/632,067; US20170306335
BI-2014/103	Broad and MIT	PCT/US2015/067138; WO2016100974
BI-2014/106	Broad and MIT	PCT/US2015/065393; WO2016094872
BI-2014/107	Broad and MIT	15/619,735; US20170349894
BI-2014/108	Broad and MIT	15/619,737; US20170349914
BI-2014/113	Broad and MIT	15/640,103; US20180112255
BI-2015/002	Broad, MIT, Harvard and Tokyo	PCT/US2016/038252; WO2016205759
BI-2015/052	Broad and MIT	PCT/US2016/038034; WO2016205613
BI-2015/053	MIT	PCT/US2016/038205; WO2016205728
10086	Broad and MIT	PCT/US2017/047458; WO2018035387
10114	Broad and MIT	PCT/US2017/053795; WO2018064208
10125	Broad and MIT	62/502,064 62/564,102
10209	Broad, Harvard, MIT, New York University and NY Genome Center	62/529,573

***Broad does not purport to grant any rights in the asterisked applicant's interest in these applications in this Agreement.**

©BEFORE OPENING OR USING THIS PRODUCT, PLEASE READ THE TERMS AND CONDITIONS SET FORTH IN THIS LIMITED USE LABEL LICENSE. YOUR USE OF THIS PRODUCT SHALL CONSTITUTE ACKNOWLEDGMENT AND ACCEPTANCE OF THESE TERMS AND CONDITIONS. If you do not agree to use this Product pursuant to the terms and conditions set out in this Limited Use Label License, please contact Promega Corporation within ten (10) calendar days of receipt to return the unused and unopened Product for a full refund; provided, however, that custom-made Products may not be returned for a refund.

Promega is selling this Product under a research use only license from Caribou Biosciences, Inc., under certain patents and patent applications owned or controlled by Caribou. For information on obtaining rights beyond those set forth in this Limited Use Label License please contact licensing@promega.com or Promega Corporation, 2800 Woods Hollow Road, Madison, WI 53711 USA, Attn: Licensing.

The purchase of this Product conveys to you, the Buyer, the NON-TRANSFERABLE right to use this Product solely for Research Use, subject to the conditions set out in this Limited Use Label License. If you wish to use this Product for any purpose other than Research Use, you must first obtain an appropriate license. This Product may not be used for any purpose other than Research Use. Your right to use this Product for Research Use is subject to the following conditions and restrictions:

1. Research Use means any use for internal research purposes, (e.g., for screening, characterization, quality control/quality assurance testing, product lot release assays, and/or toxicology purposes) for the sole benefit of you, the Buyer.
2. Research Use specifically **excludes** the following:
 - Sale or otherwise transfer of this Product, including without limitation any material that contains this Product in whole or part or any progeny, modification or derivative of this Product (collectively, "Related Materials") to any other third party; and
 - Commercial purposes and prohibited activities including use of this Product and Related Materials in (i) any sale, marketing, distribution, disposition, importation or exportation, whether for a fee or any other commercial basis; (ii) any commercial activity including, but not limited to, any use in manufacturing or bioproduction; (iii) preclinical or clinical testing or other activity directed toward the submission of data to the U.S. Food and Drug Administration, or any other regulatory agency in any country or jurisdiction where the active agent in such studies comprises this Product; (iv) use in a product for, or for any use in, human or animal therapeutics, diagnostics, or prophylactics; (v) activity in an agricultural field trial or any activity directed toward the submission of data to the U.S. Department of Agriculture or any other agriculture regulatory agency; (vi) modification of the human germline, including editing human embryos or reproductive cells; (vii) stimulation of biased inheritance of a particular gene or trait or set of genes or traits; and/or (viii) use of self-administration.
3. You, the Buyer, shall use this Product and Related Materials in compliance with all applicable laws and regulations.
4. No warranties of any kind are provided to the you, the Buyer (statutory or implied), concerning this Product or Related Materials, including without limitation, as to product quality, condition, description, merchantability, fitness for a particular purpose, noninfringement of intellectual property rights, or the absence of latent or other defects, and all such warranties are hereby expressly disclaimed.
5. Your right to use this Product and Related Materials will terminate immediately if you fail to comply with these terms and conditions.

©This Product and its use are the subject of one or more of the following patents and applications controlled by Sigma-Aldrich Co., LLC (SIGMA): Patent applications and issued patents that entered their respective National Stage from PCT International Pub. No. WO 2014/089290, including, but not limited to, the following, and substitutions, divisions, continuations, continuations-in-part, reissues, reexaminations, and extensions thereof: Australia Patent Nos. 2013355214; 2017204031; and 2018229489; Canada Patent No. 2,891,347; China Patent No. CN105142669; European Patent Nos. EP 3 138 910 B1, 3 138 911 B1, and EP 3 138 910 B1; Israel Patent No. IL238856; Singapore Patent No. 11201503824S; South Korea Patent No. 10-1844123; and U.S. Patent Application Serial Nos. 15/188,911; 15/188,924; 15/188,927; 15/188,931; and 15/456,204 (the "Patent Rights").

BEFORE OPENING OR USING THIS PRODUCT, PLEASE READ THE TERMS AND CONDITIONS SET FORTH IN THIS LICENSE AGREEMENT. YOUR USE OF THIS PRODUCT SHALL CONSTITUTE ACKNOWLEDGMENT AND ACCEPTANCE OF THESE TERMS AND CONDITIONS. If you do not agree to use this Product pursuant to the terms and conditions set out in this License Agreement, please contact Promega Corporation within ten (10) days of receipt to return the unused and unopened Product for a full refund, provided, however, that custom-made Products may not be returned for a refund.

The purchase of this Product conveys to you (the "Buyer") the NON-TRANSFERABLE right to use the Product for Licensed Research Use (see definition below) subject to the conditions set out in this License Agreement. If you wish to use this Product for any purpose other than Licensed Research Use, you must first obtain an appropriate license (see information set out below).

This Product may not be used for any purpose other than Licensed Research Use. Your right to use this Product for Licensed Research Use is subject to the following conditions and restrictions:

1. "Licensed Research Use" means any use for research purposes, with the following exceptions: (i) Buyer may not sell or otherwise transfer the Product (including without limitation any material that contains the Product in whole or part) or any Related Material to any other third party; (ii) Buyer may use the Product, components of the Product, and any Related Material, only for your internal research within the Field, and not for any Commercial Purposes (for any use of Product, components of a Product and any Related Material by any entity, including but not limited to contract service providers ("CRO") offering services for consideration, please contact Promega for additional information), and this License Agreement does not grant or in any way convey any rights to Buyer to further edit Products, components of Products or Related Materials under any claims in the foregoing patents and patent applications; (iii) Buyer shall use the Product and any Related Material in compliance with all applicable laws and regulations, including without limitation applicable human health



and animal welfare laws and regulations; (iv) SIGMA shall provide no warranties of any kind to the Buyer (statutory or implied) concerning the Patent Rights, the Product or any Related Materials, including without limitation, as to product quality, condition, description, merchantability, fitness for a particular purpose, noninfringement of intellectual property rights or the absence of latent or other defects, and all such warranties are hereby expressly disclaimed; (v) SIGMA shall expressly disclaim any warranty regarding results obtained through the use of the Product or Related Materials, including without limitation any claim of inaccurate, invalid or incomplete results; (vi) SIGMA, and their directors, officers, employees and agents, shall have no liability to the Buyer, including, without limitation, for any loss of use or profits, business interruption or any consequential, incidental, special or other indirect damages of any kind, regardless of how caused and regardless of whether an action in contract, tort, strict product liability or otherwise; (vii) the Buyer shall indemnify, defend, and hold harmless SIGMA and their current and former directors, officers, employees and agents, and their respective successors, heirs and assigns (the "Indemnitees") against any liability, damage, loss, or expense (including without limitation reasonable attorneys' fees and expenses) incurred by or imposed upon any of the Indemnitees in connection with any claims, suits, investigations, actions, demands or judgments arising out of or related to the exercise of any rights granted to the Buyer hereunder or any breach of this License Agreement by such Buyer, provided that, to the extent the foregoing is not permitted by law, the Buyer agrees, to the extent permitted by law, that it, and not the Indemnitees shall be responsible for any liability, damage, loss or expense arising out of or related to the exercise of any rights granted to the Buyer under this License Agreement or any breach of this License Agreement by Buyer, and (viii) the Product and its creation and/or use may be the subject of one or more issued patents and/or pending patent applications owned by SIGMA and the purchase of the Product and its use as defined in this License Agreement does not convey to the Buyer any license under any claims in the foregoing patents or patent applications except for the limited rights as described in this License Agreement directed toward use of the Licensed Product or commercialization thereof.

2. For purposes of Section 1 above, the following definitions shall apply:

- "Commercial Purposes" means (a) the practice, performance or provision of any method, process or service or (b) the manufacture, sale, distribution, disposition or importing of any product, in each case (a) or (b) for consideration, e.g., a fee, or on any other commercial basis.
- "Field" means use as a research tool for research purposes which includes, but is not limited to drug discovery and development, quality control/quality assurance testing and product lot release assays, provided, however, that notwithstanding the foregoing, the Field shall expressly exclude (a) any in vivo and ex vivo human or clinical use, including, without limitation, any administration into humans or any diagnostic or prognostic use, (b) the creation of transgenic rodent models and/or derivatives thereof (such derivatives including, but not limited to, rodents' cells and rodents' organs), (c) any in vivo veterinary or livestock use, or non-research agricultural use, or (d) the manufacture, distribution, importation, exportation, transportation, sale, offer for sale, marketing, promotion or other exploitation or use of, or as, a therapeutic or diagnostic for humans or animals.
- "Related Materials" means any progeny, modification or derivative of a Product.

3. Your right to use the Product will terminate immediately if you fail to comply with these terms and conditions. For information on purchasing a license (a) to this Product for purposes other than Licensed Research Use, or (b) to the foregoing patents and patent applications generally, contact your local SIGMA sales representative, who will refer you to the proper licensing representative, or in the USA call 800-325-3010.

^(e)U.S. Pat. Nos. 9,797,889, 9,797,890, 10,107,800 and 10,648,971; European patents 2970412 and 3783011; and Japanese patents 6654557 and 7280842; and other patents and patents pending.

^(f)U.S. Pat. No. 8,809,529, European Pat. No. 2635582, Japanese Pat. No. 5889910 and other patents and patents pending.

©2025 Promega Corporation. All Rights Reserved.

GloMax, Lumit, NanoBIT and NanoLuc are registered trademarks of Promega Corporation. Bio-Glo, Bio-Glo-NB and Bio-Glo-NL are trademarks of Promega Corporation.

Corning and Costar are registered trademarks of Corning, Inc. GIBCO is a registered trademark of Life Technologies Corp. GraphPad Prism is a registered trademark of GraphPad Software, Inc. JMP is a registered trademark of SAS Institute, Inc. Mr. Frosty is a registered trademark of Nalgene Nunc International Corporation. Styrofoam is a registered trademark of Dow Chemical Company.

Products may be covered by pending or issued patents or may have certain limitations. Please visit our website for more information.

All prices and specifications are subject to change without prior notice.

Product claims are subject to change. Please contact Promega Technical Services or access the Promega online catalog for the most up-to-date information on Promega products.