

**ANALYSIS OF PERCENT STUTTER FOR THE 13 CORE CODIS LOCI USING PROFILER PLUS™
AND COFILER™**

Robyn L. Ragsdale, Kelli Carter, and Petros Tsingelis

Serology/DNA Section, Florida Dept. of Law Enforcement, Tampa Regional Operations Center, Tampa, FL



The purpose of this study was to determine if the percent stutter filters used in the Kazam macro in Genotyper® were consistent with the percent stutter that we observed in our laboratory or if a modification in these values was required. Known single source samples were amplified with Profiler Plus™ and/or COfiler™ reagent kits. Amplified product was run on the ABI 310 Genetic Analyzer and analyzed using GeneScan® and Genotyper® software. Approximately 50 samples were analyzed. Percent stutter was calculated by dividing the peak height of the stutter peak (n-4) by the peak height of the called allele peak. Results for percent stutter for all alleles examined at each loci will be presented. Data will also be presented on a comparison made between the percent stutter observed for the loci found in both the Profiler Plus™ and COfiler™ kits. These results will show that our findings did not differ significantly from those in the Profiler Plus™ and COfiler™ User manuals.