

**MULTIPLEX REAL TIME FLUOROGENIC SNP ASSAYS FOR
INCREASING DISCRIMINATION OF MTDNA TESTING**

Rebecca S. Hamm, MFS; Ilona H. Letmanyi, BS; Michael D. Coble, MFS; Thomas J. Parsons, Ph.D.
Armed Forces DNA Identification Laboratory, Armed Forces Institute of Pathology, Rockville, MD



Real time fluorogenic PCR assays for allelic discrimination are widely used for single nucleotide polymorphism (SNP) testing in medical diagnostic and research settings. The TaqMan[®] allelic discrimination system that we report on here was run on an ABI 7700 sequence detection system, but is applicable to a variety of detection platforms. TaqMan[®] allelic discrimination requires amplification primers and two labeled probes for each target site. Due to the limited number of fluorescent dyes readily available, in its standard form TaqMan[®] can be “multiplexed” for only two sites. This is a limitation that serves to offset the advantages it presents in extreme sensitivity and ready availability as a well-developed commercial system.

We have devised and tested a modified approach to TaqMan[®] allelic discrimination for mtDNA typing that overcomes the multiplex limitation. This involves a multiplex amplification using the TaqMan[®] PCR primers for multiple target sites, followed by a second round of singleplex TaqMan[®] assays with the fluorogenic probes. While this is no longer a homogeneous assay, it is a simple two part assay that preserves what is by far the most important aspect of multiplexing in the forensic context: that only a single amplification needs to be performed on the casework DNA extract. We will present the design, selection strategy, and testing of 10-plex SNP assays that target sites that greatly increase the forensic discrimination of mtDNA typing when common HV1 and HV2 sequences are encountered.