

QIAGEN® MINELUTE PCR PURIFICATION INCREASES SENSITIVITY IN FOUR MULTIPLEX PCR KITS

JM Esparza, PhD and PW Wojtkiewicz, PhD

North Louisiana Criminalistics Laboratory, 1115 Brooks St., Shreveport, LA 71101

Two means determine the sensitivity of STR detection: total relative fluorescence units (RFU) and alleles detected. Several methods can be used to increase the sensitivity of low-copy number (LCN) DNA samples, such as increasing the number of PCR cycles; increasing injection time or voltage¹; increasing the amount of PCR product analyzed; and removing primers, buffers, and non-amplicon products from the PCR product. QIAGEN® MinElute PCR Purification uses a column to bind, wash and elute amplified PCR product to remove excess primers and PCR components³ leaving the sample cleaner and more concentrated. A previous study suggests that post-PCR purification of the amplified product with QIAGEN® MinElute can increase the sensitivity to effectively detect DNA profiles from <100 pg of DNA; but typical problems for LCN samples such as allele dropout, increased stutter, and sporadic contamination still exist.² However, we have shown that QIAGEN® MinElute increases the sensitivity as measured by total RFU and alleles called with four multiplex PCR kits PowerPlex® 16, AmpFISTR® SGM Plus®, AmpFISTR® YFiler™, and AmpFISTR® MiniFiler™ at the North Louisiana Criminalistics Laboratory (NLCL). A working concentration of 0.1 ng/μL of female and male DNA was used to prepare DNA samples for amplification and serially diluted two-fold to yield a final template of 1.0, 0.5, 0.25, 0.125, 0.063, 0.031, and 0.016 ng in 10 μL of TE⁻⁴ buffer. Each amount of DNA was amplified for each of the systems and subjected to QIAGEN® MinElute PCR purification. Compared to samples not cleaned QIAGEN® MinElute PCR purification, there is an overall increase in total RFU and alleles detected for DNA amounts <125 pg. However, for template amounts >125 pg, there are considerable artifacts and stochastic effects which include allele drop in, increased stutter, and pull-up from too much amplified product. Thus, the effectiveness of QIAGEN® MinElute PCR purification is dependent on the initial amount of DNA input and does not produce useable results for DNA amounts greater than 125 pg. Non-probative samples for PowerPlex® 16 and AmpFISTR® SGM Plus® amplification kits show an increase in sensitivity as measured by both total RFU and alleles detected for two sets of samples; those with no initial profile and samples with a partial profile. The use of QIAGEN® MinElute PCR purification for amplified PCR product results in an increased sensitivity for low copy number DNA samples and therefore will increase the ability to detect alleles in casework samples.