

Cleaning Log

Maxwell® 16 Instrument (Cat.# AS2000)

Serial Number:

Clean the instrument as described in the *Maxwell® 16 Instrument Operating Manual* following any spill or at regular intervals as appropriate for instrument usage. Record activities below.

Today's Date (MM/DD/YYYY): _____ **Time:** _____

User Name: _____

Please check off steps as they are completed:

- Disconnect from the Power Source.
- Remove Hardware.
- Clean Hardware.
- Clean Magnetic Rod Assembly.
- Clean Plunger Bar.
- Clean Platform.
- Clean Interior.
- Clean Exterior and Door.
- Reinsert Hardware.
- Reconnect to Power Source.

Comments: _____

Today's Date (MM/DD/YYYY): _____ **Time:** _____

User Name: _____

Please check off steps as they are completed:

- Disconnect from the Power Source.
- Remove Hardware.
- Clean Hardware.
- Clean Magnetic Rod Assembly.
- Clean Plunger Bar.
- Clean Platform.
- Clean Interior.
- Clean Exterior and Door.
- Reinsert Hardware.
- Reconnect to Power Source.

Comments: _____
