



# Cleaning Log

Maxwell® 16 Instruments (Cat.# AS3000, AS3050, AS3060)

Serial Number: \_\_\_\_\_

Clean the instrument as described in the *Maxwell® 16 Instrument Operating Manual* following any spill or at regular intervals as appropriate for instrument usage. Record activities below.

**Today's Date** (MM/DD/YYYY): \_\_\_\_\_ **Time:** \_\_\_\_\_

User Name: \_\_\_\_\_

Please check off steps as they are completed:

- Disconnect from the Power Source.
- Remove Hardware.
- Clean Hardware.
- Clean Magnetic Rod Assembly.
- Clean Plunger Bar.
- Clean Platform.
- Clean Interior.
- Clean Exterior and Door.
- Reinsert Hardware.
- Reconnect to Power Source.

Comments: \_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

**Today's Date** (MM/DD/YYYY): \_\_\_\_\_ **Time:** \_\_\_\_\_

User Name: \_\_\_\_\_

Please check off steps as they are completed:

- Disconnect from the Power Source.
- Remove Hardware.
- Clean Hardware.
- Clean Magnetic Rod Assembly.
- Clean Plunger Bar.
- Clean Platform.
- Clean Interior.
- Clean Exterior and Door.
- Reinsert Hardware.
- Reconnect to Power Source.

Comments: \_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_