

Cleaning Log

Spectrum Compact CE System (Cat.# CE1304)

Serial Number: _____

Clean the instrument as described in the Spectrum Compact CE System Operating Manual #TMD058 following any spill or at regular intervals as appropriate for instrument usage. Record activities below.

Today's Date (MM/DD/YYYY): _____

Time: _____

User Name: _____

Please check off steps as they are completed:

Comments:

Remove Consumables.

 Anode and Cathode Buffer.

 Polymer.

 Capillary.

Disconnect from the Power Source.

Remove Anode Electrode.

Clean Anode Electrode.

Clean AS Tray.

Clean Interior.

Clean Exterior and Door.

Reconnect to Power Source.

Today's Date (MM/DD/YYYY): _____

Time: _____

User Name: _____

Please check off steps as they are completed:

Comments:

Remove Consumables.

 Anode and Cathode Buffer.

 Polymer.

 Capillary.

Disconnect from the Power Source.

Remove Anode Electrode.

Clean Anode Electrode.

Clean AS Tray.

Clean Interior.

Clean Exterior and Door.

Reconnect to Power Source.

Download this form at: www.promega.com/SpectrumService